

Product datasheet for **MR220493L4V**

Camsap1 (NM_001115076) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Camsap1 (NM_001115076) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Camsap1
Synonyms:	9530003A05Rik; C77823; PRO2405
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001115076
ORF Size:	4743 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220493).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001115076.1 , NP_001108548.1
RefSeq Size:	7978 bp
RefSeq ORF:	4746 bp
Locus ID:	227634
Cytogenetics:	2 A3



[View online »](#)

Gene Summary:

Key microtubule-organizing protein that specifically binds the minus-end of non-centrosomal microtubules and regulates their dynamics and organization. Specifically recognizes growing microtubule minus-ends and stabilizes microtubules. Acts on free microtubule minus-ends that are not capped by microtubule-nucleating proteins or other factors and protects microtubule minus-ends from depolymerization. In contrast to CAMSAP2 and CAMSAP3, tracks along the growing tips of minus-end microtubules without significantly affecting the polymerization rate: binds at the very tip of the microtubules minus-end and acts as a minus-end tracking protein (-TIP) that dissociates from microtubules after allowing tubulin incorporation. Through interaction with spectrin may regulate neurite outgrowth. [UniProtKB/Swiss-Prot Function]