

Product datasheet for **MR220496L4V**

Klk11 (NM_019974) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Klk11 (NM_019974) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Klk11
Synonyms:	2310015I08Rik; hj; Pr; Prss20; T; TLSP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019974
ORF Size:	747 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220496).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019974.2 , NP_064358.1
RefSeq Size:	1285 bp
RefSeq ORF:	750 bp
Locus ID:	56538
UniProt ID:	Q9QYN3
Cytogenetics:	7 B3



[View online »](#)

Gene Summary:

This gene encodes a member of the kallikrein subfamily of serine proteases that are involved in diverse physiological functions such as skin desquamation, tooth enamel formation, seminal liquefaction, synaptic neural plasticity and brain function. The encoded preproprotein undergoes proteolytic cleavage of the activation peptide to generate the functional enzyme. This gene is located in a cluster of several related kallikrein genes on chromosome 7. [provided by RefSeq, May 2016]