

Product datasheet for **MR220505L3V**

Vpreb1 (NM_016982) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Vpreb1 (NM_016982) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Vpreb1
Synonyms:	CD179a; Vpreb-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016982
ORF Size:	429 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220505).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016982.2 , NP_058678.1
RefSeq Size:	772 bp
RefSeq ORF:	429 bp
Locus ID:	22362
UniProt ID:	P13372
Cytogenetics:	16 10.46 cM
Gene Summary:	Associates with the Ig-mu chain to form a molecular complex that is expressed on the surface of pre-B-cells. This complex presumably regulates Ig gene rearrangements in the early steps of B-cell differentiation.[UniProtKB/Swiss-Prot Function]



[View online »](#)