

Product datasheet for MR220533

B9d1 (NM_013717) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: B9d1 (NM_013717) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: B9d1
Synonyms: AW045994; B9; Eppb9
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR220533 representing NM_013717
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGCTGCAGCAAGTCCCAGCGTCTTCTGCTCATGATCACCGGGCAGGTAGAGAGCGCACAGTTCCAG
 AGTATGATGACCTCTACTGCAAGTACTGCTTTGTGTATGGCCAGGACTGGGCCCAACAGCGGGTCTGGA
 GGAGGGGATCTCACAGATAGCATCTAAGAGCCAAGATGTACGGCAAGCACTGGTGTGGAACCTCCCCATC
 GATGTGACCTTTAAAAGTACCAACCCCTATGGCTGGCCACAGATTGTACTCAGTGTCTATGGCCGGATG
 TGTTTGGGAATGATGTCGTCGAGGCTATGGAGCAGTGCATGTGCCCTCTCTCCAGGACGGCACAAAAG
 GACCATCCCCATGTTTGTGCCAGAGTCTACATCTACACTACAGAAGTTCACAAGCTGGTTTATGGGACGG
 CGACCCGAGTACACAGACCCCAAGGTGGTGGCCAGGGTGAAGGCCGGGAAGTACTCGTGTTCGCTCTC
 AGGGATTTGTTACCCTCCTCTTCAATGTGGTGACCAAGGACATGAAGAAGCTGGGCTATGATACTGGCC
 TGTGGACACACAAGGAGTCTGGGTCCTAGCCTGCCACAGGGCAACCCACAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR220533 representing NM_013717
 Red=Cloning site Green=Tags(s)

MAAASPSVFLMITGQVESAQFPEYDDL YCKYCFVYGQDWAPTAGLEEGISQIASKSQDVRQALVWNFPI
 DVTFKSTNPYGPQIVLSVYGPVDFGNDVVRGYGAVHVPLSPGRHKRTIPMFVPESTLQKFTSWFMGR
 RPEYTDPKVVAQEGREVTRVRSQGFVTL LFNVVTKDMKKLGYDTGPVDTQGVLPQGNPQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Chromatograms: https://cdn.origene.com/chromatograms/mm9043_f05.zip

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM_013717

ORF Size: 612 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_013717.3](#)

RefSeq Size: 816 bp

RefSeq ORF: 615 bp

Locus ID: 27078


UniProt ID: [Q9R1S0](#)

Cytogenetics: 11 B2

MW: 23 kDa

Gene Summary: Component of the tectonic-like complex, a complex localized at the transition zone of primary cilia and acting as a barrier that prevents diffusion of transmembrane proteins between the cilia and plasma membranes. Required for ciliogenesis and sonic hedgehog/SHH signaling. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR220533