

Product datasheet for **MR220538L3V**

Hddc3 (NM_026812) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hddc3 (NM_026812) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hddc3
Synonyms:	1110033O09Rik; C86475; MESH1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026812
ORF Size:	537 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220538).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026812.2 , NP_081088.1
RefSeq Size:	860 bp
RefSeq ORF:	540 bp
Locus ID:	68695
UniProt ID:	Q9D114
Cytogenetics:	7 D2
Gene Summary:	ppGpp hydrolyzing enzyme involved in starvation response.[UniProtKB/Swiss-Prot Function]



[View online »](#)