

## Product datasheet for **MR220539**

### **Pcp2 (NM\_001129803) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Pcp2 (NM\_001129803) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Pcp2  
**Synonyms:** L7; Pcp-2  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR220539 representing NM\_001129803  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGACCAGGAAGAAGAAAAGGGAGAGGGAGGCTCAGACCTTCTAGACAAGGCAGGTTCCACGGACCAGG  
AAGGCTTCTTCAACCTGCTGACCCACGTGCAGGGCGATCGGATGGAGGAGCAGCGCTGTTCTTGCAGGC  
TGGCCAGGCCAGAACCAGAAAGCCAGGGTGGCCCTGCTCCAGAGATGGACAATCTCATGGATATGCTG  
GTCAACACCCAGGCGCCGCATGGACGACCAGCGTGAACAGTTAATCCCTGCCTGGCTTCCAACCTA  
TCGGCCCAAGGATGGAATGCAGAAACGACCTGGGACCCTCAGCCCTCAACCCCTGCTCACCCCTCAGGA  
TCCTGCTGCACTCAGCTCCGCAGGAACAGCAGCCCCAGCCCCAGACACAAGCTCT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR220539 representing NM\_001129803  
**Red=Cloning site Green=Tags(s)**

MDQEEEKGEGSDLLDKAGSPDQEGFFNLLTHVQGDRMEEQRCSLQAGPQNPESQGGPAPEMDNLMDML  
VNTQGRMRDDQRTVNSLPGFQPIGPKDGMQKRPGLSPQPLLTPQDPAALSFRRNSSPQPQTQAP

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9050\\_d01.zip](https://cdn.origene.com/chromatograms/mm9050_d01.zip)

**Restriction Sites:** SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM\_001129803

ORF Size: 411 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001129803.2](#)

RefSeq Size: 501 bp

RefSeq ORF: 411 bp


Locus ID: 18545

Cytogenetics: 8 1.92 cM

**MW:** 14.8 kDa

**Gene Summary:** May function as a cell-type specific modulator for G protein-mediated cell signaling.  
[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR220539