

## Product datasheet for **MR220584L3V**

### Agbl2 (NM\_178755) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Agbl2 (NM_178755) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Agbl2
Synonyms:	4930524K04; A430081C19Rik; CCP2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_178755
ORF Size:	2586 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220584).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_178755.3</a>
RefSeq Size:	3462 bp
RefSeq ORF:	2589 bp
Locus ID:	271813
UniProt ID:	<a href="#">Q8CDK2</a>
Cytogenetics:	2 E1



[View online »](#)

**Gene Summary:**

Metallocarboxypeptidase that mediates deglutamylation of target proteins. Catalyzes the deglutamylation of polyglutamate side chains generated by post-translational polyglutamylation in proteins such as tubulins. Also removes gene-encoded polyglutamates from the carboxy-terminus of target proteins such as MYLK. Does not show detyrosinase or deglycylase activities from the carboxy-terminus of tubulin.[UniProtKB/Swiss-Prot Function]