

Product datasheet for MR220596

Plpp4 (NM_001080963) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Plpp4 (NM_001080963) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Plpp4
Synonyms:	C030048B12Rik; Gm1090; Ppapdc1a
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR220596 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCGGGAGCTGGCCATCGAGATCGGGGTGCGAGCCCTTCTCTTCGGGGTCTTTGTTTTACAGAAATTT
TGGATCCATTCCAGAGAGTCATCCAGCCGGAAGAGATCTGGCTCTATAAAAATCCTTTGGTGCAGTCAGA
CAACATACCTACCCGCTCATGTTTGAATTTCTTTCCTCACACCCTGGCTGTTATTTGCGTGGTGAAA
ATTATCCGGCGGACAGACAAGACTGAAATTAAGGAAGCCTTCTTAGCAGTGTCTTGGCTCTTGCTTTGA
ATGGAGTCTGCACAAACTATTAACCTAATAGTGGGAAGACCTCGCCCCGATTTCTTTTACCGCTGCTT
TCCAGATGGGGTATGAATTCGGAGATGCGCTGCACAGGTGACCTGACCTGGTGTCTGGAGGGCCGCAA
AGCTTCCCCAGCATCCATTCCTCCTTTGCCTTCTCGGGCCTCGGTTTCACAACATTCTACCTGGCTGGCA
AGCTGCACTGCTTACCAGAGAGTGGCGGGGGAAGAGCTGGAGGCTCTGTGCTGCCATCTTGCCCTTGTA
CTGTGCCATGATGATCGCCCTGTCCCGCATGTGTGACTACAAGCATATTGGCAAGATTCTTTTGTGGGA
GGAGTGATTGGCCTCATTTTGCCTATATTTGCTACAGACAACATTATCCTCCTCTGGCCAAACACAGCAT
GCCACAAACCGTATGTCAGCCTTCGAGTGCCCACTTCACTGAAGAAGGAAGAGAGGCCTACAGCCGACAG
TGCGCCAGCTTGCCCTGGAGGGGATCACCGAGGGCCCCGTA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR220596 protein sequence
Red=Cloning site Green=Tags(s)

MRELAIEIGVRALLFGVVFTEFLDPFQRVIQPEEIWLKYNPLVQSDNIPRLMFAISFLTPLAVICVVK
 IIRRTDKTEIKEAFLAVSLALALNGVCTNTIKLIVGRPRPDDFYRCFPDGMNSEMRCTGDPDLVSEGRK
 SFPSIHSSFAFSGLGFTTFYLAGKLHCFTESGRGKSWRLCAAAILPLYCAMMIALSRMCDYKHHWQDSFVG
 GVIGLIFAYICYRQHYPPLANTACHKPYVSLRVPTSLKKEERPTADSAPSLPLEGITEGPV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001080963

ORF Size: 816 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001080963.1](#), [NP_001074432.1](#)

RefSeq Size: 1088 bp

RefSeq ORF: 816 bp

Locus ID: 381925


UniProt ID: [Q0VBU9](#)

Cytogenetics: 7 F3

MW: 30.4 kDa

Gene Summary: Displays magnesium-independent phosphatidate phosphatase activity in vitro. Catalyzes the conversion of phosphatidic acid to diacylglycerol (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR220596