

Product datasheet for **MR220602L4V**

IL22ra2 (NM_178258) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	IL22ra2 (NM_178258) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	IL22ra2
Synonyms:	CRF2-10; CRF2-s1; CRF2X; IL-22bp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_178258
ORF Size:	690 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220602).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_178258.5 , NP_839989.2
RefSeq Size:	2468 bp
RefSeq ORF:	693 bp
Locus ID:	237310
UniProt ID:	Q80XF5
Cytogenetics:	10 A3
Gene Summary:	Receptor for IL22. Binds to IL22, prevents interaction with the functional IL-22R complex and blocks the activity of IL22 (in vitro). May play an important role as an IL22 antagonist in the regulation of inflammatory responses.[UniProtKB/Swiss-Prot Function]



[View online »](#)