

## Product datasheet for **MR220604**

### **Pdia4 (NM\_009787) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Pdia4 (NM_009787) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pdia4
Synonyms:	A1987846; Cai; ERp-72; Erp72
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MR220604 representing NM\_009787  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGAAGCTCCGAAAGCTTGGCTGTTGGTCTGCTCTTGGCGCTGACACAGCTGCTGGCTGCTGCGAGCG  
CCGGAGATGCACACGAAGATACTTCAGATACAGAAAATGCCACTGAGGAGGAAGAGGAGGAAGAGGAAGA  
GGATGACGATGACTTGGAAAGTTAAGGAAGAAAATGGTGTGGTCTTAAATGATGGAACTTTGATAAC  
TTTGTGGCTGACAAAGATACAGTCTACTGGAGTTCTATGCACCATGGTGTGGACTGCAAGCAGTTTG  
CTCCAGAATATGAGAAAATGCCAGTACTTTGAAGGATAATGATCCTCCCATTGCTGTAGCGAAGATCGA  
TGCAACCTCAGCATCCATGCTAGCCAGCAAATTTGATGTGAGTGGCTACCCACCATCAAGATCCTGAAG  
AAGGGACAGGCCGTTGACTATGATGGCTCCAGGACCCAGGAAGAAATTTGCCAAAGTCAGAGAAGTTT  
CCCAGCCTGATTGGACACCTCCACCTGAAGTCACTCTTTCATTGACTAAGATAACTTTGACGATGTTGT  
AAATAATGCAGACATTATTTTGGTGGAGTTCTATGCCCGTGGTGTGGACTGCAAGAACTTGCCCT  
GAGTATGAGAAAGCTGCCAAGGAGCTCAGCAAGCGCTCTCCCAATCCCTTAGCCAAAGTGGAGGCCA  
CCGAACAGACAGACCTGGCTAAGAGTTTGTATGTCTCTGGCTACCCGACTTTGAAAATTTCCGCAAGG  
AAGGCCTTTTACTACAATGGCCACGAGAGAAAATATGGGATTGTTGACTACATGATTGAGCAGTCTGGG  
CCTCCCTCAAAGGAGATCCTGACCCTGAAGCAGGTCCAGGAGTTCCTGAAGGATGGAGATGATGTGGTCA  
TCATTGGGCTCTTTCAGGGAGATGGTGACCCGGCTACTTGACAGTATCAGGATGCTGCTAACAACCTGAG  
AGAAGATTACAAGTTTACCACACTTTCCAGCCCTGAAATAGCCAAGTTCCTGAAAGTCTCCCTGGGGAAG  
TTGGTGTGACACACCCTGAGAAGTCCAGTCCAAGTATGAGCCAGGTTCCACGTGATGGATGTTTCAGG  
GCTCTACAGAAGCATCAGCCATTAAGGACTATGTAGTGAACATGCCTTGCCCTGGTGGCCATCGTAA  
GACTTCTAATGATGCTAAGCGGTACAGCAAACGCCCTGGTGGTGTATACTACAGTGTGGACTTCAGC  
TTTGATTACAGAGCTGCTACTCAGTTTTGGCGTAACAAAGTCTAGAGGTGGCCAAGGACTTCCCTGAGT  
ATACCTTCGCATTGCTGACGAGGAAGACTATGCCACGGAGGTGAAGGACCTGGGGCTCAGTGAGAGTGG  
AGAGGATGTCAATGCAGCCATCCTAGACGAGAGCGGGAAGAAAGTTTGCATGGAGCCAGAGGAGTTTGT  
TCAGACACTCTCCGGAAATTTGTCACAGCTTTTAAAAAAGGAAAATAAAGCCAGTCATCAAATCACAGC  
CAGTACCCAAGAACAACAAGGGACCAGTCAAGGTGGTGGTGGGAAAGACCTTTGATGCCATTGTGATGGA  
CCCCAAAAGGACGTCTCATTGAATCTATGCACCCTGGTGTGGGCACTGCAAGCAGCTAGAGCCTATC  
TATACCAGCTAGGCAAGAAATACAAGGGCCAGAAGGACCTGGTATCGCCAAGATGGATGCTACTGCCA  
ACGACATCACCACGACCAATATAAGGTGGAGGGCTTCCCTACCATCTATTTGCTCCCAGTGGGGACAA  
AAAGAACCAATTAATTTGAGGGTGGCAACAGAGATCTGGAGCATTTGAGCAAGTTTATAGATGAACAT  
GCCCAAAGAGGAGCAGGACCAAGGAAGAGCTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR220604 representing NM\_009787  
Red=Cloning site Green=Tags(s)

```


MKLRKAWLLVLLALQLLAAASAGDAHEDTSDTENATEEEEEEEEEEDDDLEVKKEENGVWVLDNGNFDN
FVADKDTVLLLEFYAPWCGHCKQFAPEYEKIASTLKDNDPPIAVAKIDATSASMLASKFDVSGYPTIKILK
KGQAVDYDGSRTQEEIVAKVREVSQPDWTPPEVTLSTKDNFDDVVNNADIILVEFYAPWCGHCKKLAP
EYEKAAKELSKRSPPIPLAKVDATEQTDLAKRFDVSGYPTLKIIFRKG RPFDYNGPREKYGIVDYMIEQSG
PPSKEILTLKQVQEFLLKDGDDVVIIGLFQGDGPAYLQYQDAANNLREDYKFHHTFSPEIAKFLKVS LGK
LVLT HPEKFQSKYEPRFHVMDVQGSTEASAIKDYYVKHALPLVGHRKTSNDAKRYSKRPLVVVVYSVDFS
FDYRAATQFWRNKVLEVAKDFPEYTFIAIEDEDYATEVKDLGLSESGEDVNAAILDESGKKFAMEPEEFD
SDTLREFVTAFFKKGKLPVIKSQVPKNNKGPVKVVVGKTFDAI VMDPKKDVLI EFYAPWCGHCKQLEPI
YTSLGKKYKQKDLVI AKMDATANDITNDQYKVEGFPTIYFAPSGDKKNPIKFEGGNRDLEHLSKFIDEH
ATKRSRTKEEL
    
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9005\\_b02.zip](https://cdn.origene.com/chromatograms/mm9005_b02.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_009787

**ORF Size:** 1923 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_009787.2](#), [NP\\_033917.2](#)

**RefSeq Size:** 2610 bp

**RefSeq ORF:** 1926 bp


**Locus ID:** 12304

**UniProt ID:** [P08003](#)

**Cytogenetics:** 6 22.94 cM

**MW:** 72.8 kDa

**Product images:**



Circular map for MR220604