

## Product datasheet for **MR220619L4V**

### Spag11a (NM\_153115) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Spag11a (NM_153115) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Spag11a
Synonyms:	9230111C08Rik; Bin1b; EP2e; EP2Q; Spag11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_153115
ORF Size:	207 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220619).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_153115.1</a> , <a href="#">NP_694755.1</a>
RefSeq Size:	435 bp
RefSeq ORF:	210 bp
Locus ID:	78128
UniProt ID:	<a href="#">Q8K4N2</a>
Cytogenetics:	8 A1.3
Gene Summary:	May have antimicrobial activity. May also play a role sperm maturation, storage, and protection.[UniProtKB/Swiss-Prot Function]



[View online »](#)