

Product datasheet for **MR220622L3V**

Hint2 (NM_026871) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hint2 (NM_026871) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hint2
Synonyms:	1190005L05Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026871
ORF Size:	489 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220622).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026871.1 , NP_081147.1
RefSeq Size:	636 bp
RefSeq ORF:	492 bp
Locus ID:	68917
UniProt ID:	Q9D0S9
Cytogenetics:	4 A5
Gene Summary:	Hydrolase probably involved in steroid biosynthesis. May play a role in apoptosis. Has adenosine phosphoramidase activity (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)