

Product datasheet for **MR220672**

Wbp11 (NM_021714) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Wbp11 (NM_021714) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Wbp11
Synonyms:	2510026P17Rik; D6Wsu113e; Npwbp; SIPP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR220672 representing NM_021714
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGACGAAGATCAACATCCTCCACCAAGAGTGGGAAGTTTATGAACCCACAGACCAAGCCGAAAGG
AAGCCCGAAAAAGAGAAGTGA AAAAGAACAAAAACAGCGCATGATGGTCCGGGCTGCAGTTTTGAAGAT
GAAAGACCCCAAACAAATTATCCGGGACATGGAGAAATTGGATGAAATGGAGTTTAAACCCAGTGCAGCAG
CCACAGTTGAATGAGAAAGTACTGAAAGACAAGCGAAAAAAGCTGCGGGAAACCTTTGAACGCATTCTAC
GCCTTTATGAAAAAGAAAACCCAGATATTTACAAAGAAGTGAAGAAAGTGGAAAGTGAATATGAACAGAA
GAGAGCTCAGCTGAGCCAGTACTTTGATGCTGTTAAGAATGCTCAGCATGTGGAAGTAGAGAGCATCCCT
TTGCCAGACATGCCGCATGCTCCTTCCAACATCTTGATCCAGGACATTCTCTTCTCGTGCACAGCCCC
CCTCCATCCTGAAGAAAACCTCAGCGTATGGGCCCTCAGCTCGGGCCGTGTCTATCTTCTCTCCTGGG
ACATGGTGTTCCTCGCTTGCCTTGGCAGAAAACCGCCCGCCCTCCACCTGGCCCCCTCCTCCTCAG
GTTCTGCAGATGTACGGCCGAAAGTGGGCTTCGCTCTAGATCTTCCCTTCGGAGGCGAGATGAGGACA
TGCTGTATAGTCCAGAGCTTGCTCAGCGGGTCAATGACGATGACATGTCCAGCACCAGTGAAGACGACGG
CTATCCTGAGGACATGGATCAGGATAAGCATGATGACAGTACTGAGGACAGTGCAGTGCAGAGTGCAGT
GCAGAGAGTGACGGGATGAGTTCGGGCACCGCAGGACAGTGCAGCGCAACACTGAGGAGAAGAAGT
CAGGTCTAAGCGTCAGATTTGCTGATATGCCTGGGAAATCAAGGAAGAAGAAGAAAAACATGAAGGAGCT
GACTCCTCTTCAAGCCATGATGCTGCGAATGGCAGGTGAGGAAATCCCTGAGGAAGGACGGGAAGTGGAG
GAATTTTCAAGAGGAGGAGGATGCCGACGACTCGGATGATTCTGAAGCCGAGAAGCAGTGCAGAAAGCAGC
ATAAAGACGATGGCCACTCTGACAGCACGGTGCAGCCTTCCCAGCAGCAGGCTCCTCCACAGTCTGC
TCCTGCTTCTCAGATACAAGCCCTCCCATGCCAGGACCACCACCTCTTGACCACCACCTGCCCCACCC
TTACGGCCTCCTGGACCCTACAGGCCTTCTCCTGGGCCCTCCAGGGCTCCTCCATTCTGAGAC
CACCTGGAATGCCAGGAATCCGAGGGCCTTACCACGACTTTTACCTCCAGGACCACCAGGCGGACC
TCCTGGCCCTCCCCAGGCCACCTCCAGGTCTGCCCCCTGGACCCCTCCTCGGGGACCCCAACCAAGG
CTACCTCCTCCTGCACCTCCAGGTATCCCTCCTCCTGCTCCTGGCATGATGCGCCACCCTTGGTTCCTC
CTTTGGACCTGCCCCGCTGGCCTTTTCCACCAGCTCCCTTGCCCAACCCAGGGGTGTTGAGTGGCCC
ACCCAGCCTCATTGAGCGACCGAAGGCGGATGATGCAAGTGCCGCCACCATCGAGAAGAAAGCCACAGCG
ACCATCAGTGCCAAAGCCAGATCACCAATCCCAAGGCAGAGGTACACGTTTTGTGCCACTGCACTGA
GGGTACGCCGAGAGAATAAGGGGCTACTGCTGTTCCCAAGAAAGGTGAGAGGATGACTCTGCTGTGCC
TGTTGCCAAAGCAGCCCCGAGATCTGGTCTTCAAGTGTGCTCAGTACAAACTAAGGACGATGTTTAT
GAAGCCTTTATGAAAGAGATGGAAGGGCTGCTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR220672 representing NM_021714
 Red=Cloning site Green=Tags(s)


```
MGRRSTSTKSGKFMNPTDQARKEARKRELKKNKKQRMVRAAVLKMKDPKQIIRDMEKLDMEFNPVQQ
PQLNEKVLKDKRKLRETFERILRLYEKENPDIYKELRKLEVEYEKRAQLSQYFQAVKNAQHVESIP
LPDMPHAPSNI LIQDIPLGAQPPSILKKT SAYGPPARAVSILPLLGHGVPRLPPGRKPPGPPGPPPPQ
VLQMYGRKVGFALDLPPRRRDEMLYSPELAQRGHDDMSSTSEDDGYPEDMDQDKHDDSTEDSDTDRSD
AESDGDEFHREDSERDNTEKKSGLSVRFADMPGKSRKKNMKEKLTPLQAMMLRMAGQEIP EEGREVE
EFSEEDADSDSSEAEKQSQKHKDDGHSSTAAASSQQAPPQSAPASQIQAPPMPGPPPLGPPPAPP
LRPPGPPTGLPPGPPGAPPFLRPPGMPGIRGPLRLLPPGPPGRPPGPPGPPPLPPGPPPRGPPPR
LPPAPPGI PPPRGM RPPLVPLGPAPPGLFPPAPLPNGVL SAPPSLIQRPKADDASAATIEKKATA
TISAKPQITNPKAEVTRFVPTALRVRRENKGATAVPQRSEDDSAVPVAKAAPRSGPSVAVSVQTKDDVY
EAFMKEMEGLL
```

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9089_h01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_021714

ORF Size: 1923 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_021714.4](#), [NP_068360.4](#)

RefSeq Size: 2722 bp

RefSeq ORF: 1926 bp

Locus ID: 60321


UniProt ID: [Q923D5](#)

Cytogenetics: 6 66.72 cM

MW: 69.9 kDa

Gene Summary: Activates pre-mRNA splicing. May inhibit PP1 phosphatase activity.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR220672