

Product datasheet for **MR220685L3V**

Tmprss3 (NM_080727) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmprss3 (NM_080727) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tmprss3
Synonyms:	MGC130589
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_080727
ORF Size:	1359 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220685).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080727.2 , NP_542765.2
RefSeq Size:	2738 bp
RefSeq ORF:	1362 bp
Locus ID:	140765
UniProt ID:	Q8K1T0
Cytogenetics:	17 A3.3
Gene Summary:	Probable serine protease that plays a role in hearing. Acts as a permissive factor for cochlear hair cell survival and activation at the onset of hearing and is required for saccular hair cell survival. Activates ENaC (in vitro).[UniProtKB/Swiss-Prot Function]



[View online »](#)