

Product datasheet for **MR220686L4V**

Rnd2 (NM_009708) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rnd2 (NM_009708) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rnd2
Synonyms:	Arhn; rho7; Rohn
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009708
ORF Size:	681 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220686).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009708.1
RefSeq Size:	684 bp
RefSeq ORF:	684 bp
Locus ID:	11858
UniProt ID:	Q9QYM5
Cytogenetics:	11 D
Gene Summary:	May be specifically involved in neuronal and hepatic functions. Is a C3 toxin-insensitive member of the Rho subfamily (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)