

Product datasheet for **MR220708L3V**

Cdk11b (NM_007661) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cdk11b (NM_007661) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cdk11b
Synonyms:	AA989746; Cdc2l1; Cdc2l2; Cdc11b; Cdk11; CDK11-p46; CDK11-p58; CDK11-p110; p58
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007661
ORF Size:	2352 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220708).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007661.3
RefSeq Size:	2589 bp
RefSeq ORF:	2355 bp
Locus ID:	12537
UniProt ID:	P24788
Cytogenetics:	4 86.73 cM
Gene Summary:	Plays multiple roles in cell cycle progression, cytokinesis and apoptosis. Involved in pre-mRNA splicing in a kinase activity-dependent manner. May act as a negative regulator of normal cell cycle progression.[UniProtKB/Swiss-Prot Function]



[View online »](#)