

## Product datasheet for **MR220745**

### **Sppl2c (NM\_199019) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Sppl2c (NM_199019) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Sppl2c
Synonyms:	4933407P14Rik; IMP5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MR220745 representing NM\_199019  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGCATGCCTGGGTTCCCTCCACCCCTTGGGCTCACTCCTCCTCTTCTCCTCCTCCTCCTCAGCC  
 CCGAGGCCCGGGGAGAATACGGCCTGGTCCGAGTCGTGTCAAAGAAGTGGAGCAAGGACTACTGTGTCCT  
 GTACAGCTCTGACTACGTCAACCTCCCCGGGACCTGCACCATGCCCCACTTCTGTCCCTGCATGATGGC  
 ACCAAGACACCCTGGTGTCCAGATGAAGACTCTTCCACCAGGCCAGGACAGTCCCTAGACAGCGAC  
 CCCTGCACCAGACTACTACCATGGTCACGCGTGGCAACTGTAGCTTCTATGCCAAGGGCTGGCTGGCTCA  
 AGATCAAGGGCCCGAGGGTTGCTCATTGTGAGCCGGCCAGGAATCAACAGTGTTCAGACACCATCTCA  
 AAGCCCCAGGACCCAGCAAGCCCTGGCCTGCCCTCACCATCCCAGTGGCCGTGCTCCGCTACACTGACA  
 TGTTGGACATTGTCAGCCACACTTACGGTGACACCGATGTCGAGTAGCCATGTTTGACCCCTCGAGCC  
 TGTTACTGACTATAACATGGCGATTATCTTCATCTTGGCAGTGGGTACTGTTGCTGCAGGTGGCTACTGG  
 GCTGGTTTTGATGGAGGCAACAAGCTGCAGAGGCCAGGCCAAAGAGGAGGGGGCTTGGTGGTCACA  
 ATCAACAACAACAGTGGCAGCTGAGAGGTCCCAAAGGGCCTGGGAAGATGACGACTTTGAAGATGCACC  
 TATGGATTTACACCAGCCATGACGGGCGCTGTGGTCACCATGTCCTGCTCCATCATGATCCTGCTGAT  
 TTTTTCTATGACTGCTTTGTCTACGTATGATTGGGATCTTAGCCTGGGTGCCAGTACCGGCCCTTACA  
 GCTGCCTGGCACCACTACTGTGCCACCTGCCCTATGGGATACCAAGTGGTCTTGCCTGGACAACGGGT  
 TTCTGTAACATGGCCATTACTACTGTGGCAGGCTGTGTGCCATGGTACTGCTCCTCTGGGTCACTCCAC  
 CGAAATGAGGACCATGGGCATGGCTCCTGCAGGACACTTTGGGTGTCGCTACTGCCTTTTGTCTGTA  
 GCGGGTACGGCTGCCACGTTCAAGAAGTGCCTTTGTTCTGCTGGCCCTACTAGCCTTCGATGTCTT  
 CTTTGTCTTTATCACTCCCTGTTCAACAAGACTGGTGAGAGCATTATGGTAGAAGTGGCATCAGGTCCA  
 GCAGATTCTTCAAGCCACGAGAGGTTACCAATGGTGTCAAGGTGCCCGGCTAAGCTTTTCCGCCTTGA  
 CGTTGTGCAACCAGCCCTTCTCCATCCTTGGCTTTGGGACATTGTTGTCCTGGCTTCTAGTGGCCTA  
 CTGTCACCGCTTCGACATGCAAGTCCAATCTCGCCAGGTCTACTATATGGCCTGTACAGTGGCCTATGCC  
 GTGGGCTGCTGGTACCTTTGTCGCCATGATCCTCATGCAGATGGGCCAGCCTGCCCTGCTCTACTTAG  
 TGTCCAGCACTCTGCTCACCAGCCTTGTGTGGCCACCTGCCCAAGAGTTCACCCTTTTCTGGACTGG  
 CCAGGGCAGAGCCAAGATTCTGCTGAGCCTGTGGCACAGCCCTGCATTGCCTGCTGTGGGCTCCAAG  
 ATGAAGCTGGAGGATGCGAAGGATTCCTGTACAACCAACAGGTTTGGAGCAGGCCGTTGATGGAGAGTCAG  
 GGGACTTAGAGAGCAGCACTGGGACGATATGGCAGAGATGGTACCCTATCTGAGGATGAAGCCACGAG  
 TCCAGAGGGCCACAGCGAGATTCTGAGGGCTGGAGTGACACCAACCTGGACCCCAATGAGCTACCTCT  
 GGCTCCCCATGGCCCTAGAGGCCATGTTGATCCCACTGATCCAACCAATCCCACACCCCTCCGAGCTGG  
 GCCACATCCGCACTCAATCTCGGGTCCATGATAGCAGCCTGCCCTGGATGGGGTTGCACAAGAGAAAGGG  
 CTTGAAGGTGAAGAAGAGTATGTCAGCCAGGCTCCCCTA

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
 TGGATTACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR220745 representing NM\_199019  
 Red=Cloning site Green=Tags(s)

```


MACLGSLHPLGSLLLLFLLLLLSPEARGEYGLVRVVSKNWSKDYCVLYSSDYVNLPRDLHHPALLSLHDG
TKTPWCPDEDSFHQAQDSSPRQRPLHQTTMTVTRGNCSFYAKGWLAQDQGAQGLLIVSRARNQQCSDTIS
KPQDPSKPWPAL TIPVAVLRYTDMLDIVSHTYGDTDVVRVAMFAPLEPVTDYNMAIIFILAVGTVAAGGYW
AGLMEANKLQRRQAQRGGGLGGHNQQQTVAASRSQRAWEDDDFEDAPMDFTPAMTGAVVTMSCSIMILLY
FFYDCFVYVMIGIFSLGASTGLYSCLAPILCHLPLWRYQWVLPQQRVSVTWPLLLL LAGLCAMVTVLWVIH
RNEDHWAWLLQDTLGVAYCLFVLRVRVRLPTFKNCTFLLLALLAFDVFFVFITPLFKTGESIMVEVASGP
ADSSSHERLPMVLKVPRLSFSAL TLCNQPFILGF GDIVVPGFLVAYCHRFDQVQSRQVYVYMACTVAYA
VGLLVTFVAMILMQMQPALLYLVSSTLLTSLAVATCRQEFTLFTWGQGRAKIPAEPVAQPCIASAVGSK
MKLEDAKDSRTTNRFEQAVDGEDSGDLESSTGDDMAEMVTLSEDEATSPEGHSESSEGWSDTNLDPNELPS
GSPMALEAMLIPLIQPIPHSELGHIRTQSRVHDSLPLWMLHHRKGLKVKKSMSAQAPL
  
```

SGP TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9004\\_g05.zip](https://cdn.origene.com/chromatograms/mm9004_g05.zip)

**Restriction Sites:** SgfI-RsrII

**Cloning Scheme:**




**ACCN:** NM\_199019

**ORF Size:** 2070 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u>NM_199019.2, NP_950184.2</u>
<b>RefSeq Size:</b>	2225 bp
<b>RefSeq ORF:</b>	2073 bp
<b>Locus ID:</b>	237958
<b>UniProt ID:</b>	<u>A2A6C4</u>
<b>Cytogenetics:</b>	11 E1
<b>MW:</b>	76.7 kDa
<b>Gene Summary:</b>	Intramembrane-cleaving aspartic protease (I-CLiP) that may be able to cleave type II membrane signal peptides in the hydrophobic plane of the membrane.[UniProtKB/Swiss-Prot Function]

**Product images:**


Circular map for MR220745