

Product datasheet for **MR220776L3V**

Dcst1 (NM_029974) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dcst1 (NM_029974) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dcst1
Synonyms:	A330106H01Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029974
ORF Size:	2196 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220776).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029974.1 , NP_084250.1
RefSeq Size:	2312 bp
RefSeq ORF:	2199 bp
Locus ID:	77772
UniProt ID:	Q059Y8
Cytogenetics:	3 F1
Gene Summary:	This gene encodes a transmembrane protein. A similar gene in human has been associated with dendritic cells. [provided by RefSeq, May 2015]



[View online »](#)