

Product datasheet for **MR220784**

Tnfaip2 (NM_009396) Mouse Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Tnfaip2 (NM_009396) Mouse Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | Tnfaip2 |
| Synonyms: | B94; Exoc3l3; M-sec; tnfb94; Tnfip2 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

**ORF Nucleotide
Sequence:**

>MR220784 representing NM_009396
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGCTGAAGATGGTGACTTTCTTCCAAGGGTTCTCGCCAGCAGTCTGTGCCAGGGACCCTGAACCTCG
 CTGTAAGCCCCCAAAGCCACGCTCTACATCCGAGGCTGAATCAGAGACCTCCATGTCTGAGGCGTCCTC
 TGAGGACCTGATGCCATCCCGGAGGCTCCCGATGGGGAGGAGGAGTCTGCGAAGAAGAAGGAGAAGAAG
 TCCAAAGGACTGGCCAACATGTTCAAGTGTCTTACCAAAGGGAAAAAGAAAAAAGGACCAGCCAGAT
 TATCAGATCTGGAAGTGCAGCCCAAGCCAGGCCGAGTTAGATGGTCCACTGCCACAGTGAGGAGCT
 CAAGGAGGCCCTGGAGCACGGCGGCTGGAGGTGGCTGGCAGGTGCTGGCGTTGGAGAGGCAGCTCGAG
 GCTGCGGCAGCGCGGGTGGCATGAGCAACGAGGAGCTGGTGTGGCGGCAGAGCAAGGTGGAGGCTCTGT
 ACGTGCTGCTTTGCGACCAGGTGCTCGGGTACTGCGGCGGCCGCTGGAGGCGGCCCGAGCGGCTGAG
 CCAGGCGCTAGCCGTGGTGTGCGAGGAGGAGCTCGAGGACCGCGGGCGTCCGGGGGACCCTTGGCGGCG
 GCACTGGAGGCCACGCGCCCGCGGCGATGGCTGCAGCGGTGGAGGGGCGTCTGGCGGAGGTGGCAGCGG
 AACGTCTGGACGCGCAGCCGGCCACAGCGCCCGAGGGCCGCTCGAGGCGGAGAGCAGGTTCTGCACAT
 GGGTCGCACCATGAAGGAGGACCTGGAGGTGGTGGTGGAGAGGCTGAAGCCGCTGTTCCCGACGAGTTC
 AACGTCGTGCGCACCTACGCCGAAAGCTACCACTACCACTTCGCCAGCCACCTGTGCGCCTTGGCGCAGT
 TCGAGTTGTGTGAAAGGGATACCTACTTGCTGCTGCTCTGGGTGCAGAACCTCTACCCCAATGATATTCT
 CAACAGCCCGAAGTTGGCACAGGAGCTTCAAGGTGTGGGCTTGGGAGCCTCTGCCCCAAAGCAGATC
 AGATTGCTGGAGGCCATGTTCTTGTCCAATGAGGTGACCAGCGTGAAGCAGCTCATGGCCCGAGCTTTAG
 AGCTGGAGTCTCAGCGCTGGACCCAGGATGTGGCTCCCCAGAGCCTGGATGGCCACTGTACAGTGAGCT
 GGCCATCGACATTCTCCAGATCATTTACAAGGCCAGACTAAGGCTGAGAACATCACTTCGGACGTGGGG
 ATGCAGATAAAGCAGTTGCTGCTGGTGGAGCTGGCTGCACTACTGAGGAGCTACCAGCGCCCTTTGATG
 AGTTTCTAGAGAAGAGCAAACTCCTGAGGAATTACAGAGTCAACATCATGGCCAACATCAACAAGTCCCT
 GTTCTTCTGGACTTCCGTGGAGCAGAAATGGCAGATATCTCATGATTCTCTGAACCGCTGCTGGAGCCT
 TTGAAAGACCTCAAGGCCATGGCTTCGACACCCTGCTCCAGAGCCTATTTTGGACCTGAAGCCGCTGT
 TTAAGAAGTTCACACAGACCCGTTGGGCAAATCCAGTTGAGACCCTGGAGGAAATCATCACTACTGTGAG
 CAGCAGCCTGCCTGAATTCTCAGAGCTGCAGGACTGTTTCCGGGAGGAGCTCATGGAGACTGTGCACCTG
 CACCTAGTGAAAGAATACATCATCAGGCTCTGCAAACGGCGCCTGGTCTCAAGACGGCCGAGCAGCAGC
 AGCAGCTGGCAAGGCACATCCTGGCCAACGCGGATGCCATTACAGGGTTCTGCACTGAGAATGGCTCCAC
 TGCAACCTGGCTGCACCGCGCCCTCCCTATGATTGCTGAGATCATTGCGCTGCAAGATTCAGTGCCATT
 AAGATCGAGGTGGCCACATACGCCACTTGGTACCTGACTTCAGCAAAGGCCACCTGAACGCCATCCTGG
 CCATCAAGGGCAATCTACCAAGCAGTGAAGTCAGGAGCATCCGGAACATACTGGACATCAACACGGGAGT
 GCAGGAGCCTCCAGGCCCTGTTCTCCCTTATAAAGTTACT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR220784 representing NM_009396
Red=Cloning site Green=Tags(s)


MLKMYTTFQGFPGQSSVPGTLNFAVSPQKPRSTSEAESETSMSEASSEDLMPSPEAPDGEESAKKKEKK
SKGLANMFSVFTGKKKKKDKPRLSDLEVQPKPRPELDGPLTVEELKEALEHGRLEVAWQVLALERQLE
AAAAAGGMSNEELVWRQSKVEALVYLLCDQVLGVLRRPLEAAPERLSQALAVVSQEELDRRASGGPLAA
ALEATRRRWLQRWRGVVAEVAERLDAQPATAPTEGRSEAESRFLHMGRTMKEDLEVVVVERLKPLPDEF
NVVRTYAESYHYHFASHLCALAQFELCERDTYLLLLWVQNLYPNDILNSPKLAQELQGVGLGSLPPKQI
RLLEAMFLSNEVTSVKQLMARALELESQRWTQDVAPQSLDGHCHSELAIDILQIISQGQTKAENITSDVG
MQIKQLLLVELAALLRSYQRAFDEFLEKSKLLRNYRVNIMANINCLFFWTSVEQKWQISHDSLNRLLLEP
LKDLKAHGFDTLLQSLFLDLKPLFKKFTQTRWANPVETLEEIITTVSSSLPEFSELQDCFREELMETVHL
HLVKEYIIRLCKRRLVLKTAEQQQQLARHILANADAIQGFCTENGSTATWHLRALPMIAEIIRLQDSSAI
KIEVATYATWYPDFSKGHLNATLAIKGNLPSSEVRSIRNILDINTGVOEPPRPLFSLIKVT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9101_g04.zip

Restriction Sites: Sgfl-Mlul

Cloning Scheme:




ACCN: NM_009396

ORF Size: 2073 bp

| | |
|-------------------------------|--|
| OTI Disclaimer: | <p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p> |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_009396.2 , NP_033422.2 |
| RefSeq Size: | 3590 bp |
| RefSeq ORF: | 2076 bp |
| Locus ID: | 21928 |
| UniProt ID: | Q61333 |
| Cytogenetics: | 12 60.99 cM |
| MW: | 78.1 kDa |
| Gene Summary: | May play a role as a mediator of inflammation and angiogenesis.[UniProtKB/Swiss-Prot Function] |

Product images:



Circular map for MR220784