

Product datasheet for **MR220814L3V**

PPP2R3A (NM_172144) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PPP2R3A (NM_172144) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	PPP2R3A
Synonyms:	3222402P14Rik; A730042E07
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_172144
ORF Size:	1590 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220814).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172144.3 , NP_742156.2
RefSeq Size:	4778 bp
RefSeq ORF:	1593 bp
Locus ID:	235542
Cytogenetics:	9 E4



[View online »](#)