

## Product datasheet for **MR220824L3V**

### Phb (NM\_008831) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Phb (NM_008831) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Phb
Synonyms:	Bap32
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008831
ORF Size:	816 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220824).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008831.4</a> , <a href="#">NP_032857.1</a>
RefSeq Size:	1815 bp
RefSeq ORF:	819 bp
Locus ID:	18673
UniProt ID:	<a href="#">P67778</a>
Cytogenetics:	11 D
Gene Summary:	Prohibitin inhibits DNA synthesis. It has a role in regulating proliferation. As yet it is unclear if the protein or the mRNA exhibits this effect. May play a role in regulating mitochondrial respiration activity and in aging (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)