

Product datasheet for MR220840

Nkain2 (NM_001025286) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nkain2 (NM_001025286) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Nkain2
Synonyms: 6330571D19Rik; AW455467; Tcba1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR220840 representing NM_001025286
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGGTTATTGCAGTGGCAGGTGCACGCTTATCTTTATCTGTGGCATGCAACTGGTTTGTGTGCTAGAGA
 GGCAAATATTTGATTTCTTGGATATCAGTGGCACCTATCCTGGCAAATTTGTACATATAATCATCGT
 CATCCTTGGGTTGTTTGGAACTATTCAATATAGACCTCGTTACGTAACAGGATATGCTGTCTGGCTAGTT
 CTCTGGGTGACCTGGAATGTATTTGTTATCTGCTTCTATTTGGAGGCTGGAGATCTCTCAAAGGAAACAG
 ACCTTATCCTGACATTTAATATTCAATGCACCGATCTTGGTGGATGGAGAATGGACCAGGATGTATGGT
 GACGTCTGTGACACCTGCCCCAGACTGGGCTCCAGAAGACCATCGCTACATCACTGTCTCAGGGTGT
 CTGGACTACCAATACATAGAAGTGGCTCATAGTTCCCTCCAGATTGTCCTCGCAGTAAGTGTGACAACA
 AACCACTCCTTCTAGTGTCTGGAATGTTTTTCTGGGCAAACACTGTGTGCAAAAGGACTGTATGTT
 GATGTGAATGTGTGTTGGCATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR220840 representing NM_001025286
 Red=Cloning site Green=Tags(s)

MGYCSGRCTLIFICGMQLVCLERQIFDFLGYQWAPILANFVHIIIVILGLFGTIQYRPRVYVTGYAVWL
 LWVTWNVFVICFYLEAGDLSKETDLILTFNISMHRSWWMENGPVMTSVTPAPDWAPEDHRYITVSGCF
 LDYQYIEVAHSSLQIVLAVSVDNKPLLLVFLFLFGLGKHCYQKDCMLMWNVCGHFF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001025286

ORF Size: 588 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001025286.2](#), [NP_001020457.2](#)

RefSeq Size: 1233 bp

RefSeq ORF: 591 bp

Locus ID: 432450

UniProt ID: [Q4PNJ2](#)

Cytogenetics: 10 A4

MW: 23.1 kDa