

Product datasheet for **MR220850L3V**

0610010K14Rik (NM_001177607) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	0610010K14Rik (NM_001177607) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	0610010K14Rik
Synonyms:	1110020A23Rik; AL033328; AU045833; Bap18
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001177607
ORF Size:	327 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220850).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001177607.1 , NP_001171078.1
RefSeq Size:	725 bp
RefSeq ORF:	330 bp
Locus ID:	104457
UniProt ID:	Q9DCT6
Cytogenetics:	11 B3
Gene Summary:	Component of chromatin complexes such as the MLL1/MLL and NURF complexes. [UniProtKB/Swiss-Prot Function]



[View online »](#)