

## Product datasheet for **MR220863L3V**

### **Nktr (NM\_010918) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Nktr (NM_010918) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nktr
Synonyms:	5330401F18Rik; D9Wsu172e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010918
ORF Size:	4359 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220863).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_010918.2</a> , <a href="#">NP_035048.3</a>
RefSeq Size:	7026 bp
RefSeq ORF:	4362 bp
Locus ID:	18087
UniProt ID:	<a href="#">P30415</a>
Cytogenetics:	9 72.57 cM
Gene Summary:	PPIase that catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and may therefore assist protein folding. Component of a putative tumor-recognition complex involved in the function of NK cells.[UniProtKB/Swiss-Prot Function]



[View online »](#)