

Product datasheet for MR220875

Olfr1010 (NM_207149) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr1010 (NM_207149) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr1010
Synonyms:	MOR213-1; MOR213-1_p; MOR213-7P; Olfr1011; Olfr1011-ps1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR220875 representing NM_207149 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGATGAGGAAAATATTACTGAAGTGACAGAGTTTATTCTTTTAGGATTCACAGAAGATGCGGTGTTAC
AGCAAGTACTCTTTTCATCTTCCTTACCATTTACATCATTAGTCTGCTAGGGAATTTACTCTGATTTTC
TCTCATCTGTGCTGATTCAGACTCCATACCCCATGATTTCTTCATTGGAAACCTGTCATTCTGGAT
CTGTGGTATTCTCTGTCTATGCTCCAAAATCCTGATAACATGCATCACTGAAGACAAAAGCATCTCCT
TTGCTGGTTGCCTGGCTCAGTTCTTCTCTGCTGGGCTGGCCTATAGTGAATGCTACCTACTGGCTGC
CATGGCTTATGATCGTTATGTGGCCATCTCCAACCCCTTGCTTTACTCCCAGGCTATGCTCCAAGGTTA
CGTGCCAGTCTTGTTCGAGCATCTATCTTGGTGGCTTTCGAAACTCCACCATTATCACCAAGTACAT
TTACCCTGAGTTTCTGTGAGACAATATCATTGATGACTTCTTCTGTGATCTGCCTCCTCTGTGAAGTT
GGCATGTGATGTGAAAGAGAGCTACCAGGCTGTGCTCTATTTTCATATTGGCTTCTAATGTCACTCCC
ACTGTGCTATTTTGGCCTCCTACCTCTTCATCATTGTAGCCATTCTGAAGATCCGCTCCACCCAGGGCC
GCCTCAAGGCCTTCTCTACATGTGGTCCCACCTGACAGCTGTCACTTTGTAATGTTCAATCCTCTT
TATTTACTCCCGCCAAGCACAAGTTATGCCCTGGAGAGGGACAAAGTGGTGTCACTGTTCTATACTGTA
GTGATTCGAATGCTTAACCCATATGATCTACAGCTTACGAAATAAAGATGTCAAGGATGCATTGAGGAAAA
TGTTAGATACAGTCAAGTTGTCA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR220875 representing NM_207149
Red=Cloning site Green=Tags(s)

MDEENITEVTEFILLGFTEDAVLQQVLFIFLTIYIISLLGNITLISLICADSRSLHTPMYFFIGNLSFLD
 LWYSSVYAPKILITCITEDKSISFAGCLAQFFF SAGLAYSECYLLAAMAYDRYVAISNPLLYSQAMSPRL
 RASLVAASYLGGFANSTIITSETFTLSFCRDNIIDFFCDLPLVLKACDVKESYQAVLYFILASNVITP
 TVLILASLYLFIIIVAILKIRSTQGRKAFSTCGSHLTAVTLYYGSILFIYSRPSTSYALERDKVVSVFYTV
 VIPMLNPMIYSLRNKDVKDALRKM LDTVKLS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9083_d08.zip

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM_207149


ORF Size: 933 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_207149.2 , NP_997032.2
RefSeq Size:	936 bp
RefSeq ORF:	936 bp
Locus ID:	258255
Cytogenetics:	2 D
MW:	34.9 kDa
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR220875