

Product datasheet for **MR220908L4V**

Atp5g3 (NM_175015) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Atp5g3 (NM_175015) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Atp5g3
Synonyms:	6030447M23; Atp5mc3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_175015
ORF Size:	426 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220908).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_175015.2 , NP_778180.1
RefSeq Size:	729 bp
RefSeq ORF:	426 bp
Locus ID:	228033
UniProt ID:	P56384
Cytogenetics:	2 C3


[View online »](#)

Gene Summary:

The protein encoded by this gene is a subunit of mitochondrial membrane ATP synthase, the enzyme that catalyzes ATP synthesis during oxidative phosphorylation. This gene encodes subunit 9, which is present in multiple copies in the transmembrane part of the ATP synthase complex. Phenotype and gene expression profiles suggest correlations between this gene and alcoholism- and obesity-related phenotypes. Alternative splicing results in multiple transcript variants and protein isoforms. [provided by RefSeq, Sep 2014]