

Product datasheet for **MR220920**

Ccdc169 (NM_177203) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ccdc169 (NM_177203) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ccdc169
Synonyms: A730037C10Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR220920 representing NM_177203
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGAGAGGGCCACGGAGACACTTTCGAGGGGTGAGCACAGACCGCCTGAAGCTGGAGCTGCTGGAGG
AAATCCACATGAAGGATGTTGTGCAGCTCTCCACACTTGAGATACGGCACAAGATAGCAGAGCTGGAGGC
CAACCTCAACGGTGACCTTGCAGGCAGTGAGTGGAAAACCGGTACGAGACCCAGCTGGAGCTGAATGAC
CAATTGAAAAGCAAATAGTGTCTCTGAAAGAGAAAATGGAAAAGATGCGTGGAAACCTTCAGATAGAC
TGTCTTCCATTCTGTCTATGAGAAGATGCCAGTGGAACTTTAAACGTATTACTTAAACAAGTAAAAA
AGAAAAAGGAGTCTCGAAAGCCAAGTGAAGAATATGCGTTTAGACTGGAGCAGGAGTCAAAGGCATAC
CACAGGACCAACAATGAGCGCCGTTCTACATCGCGGAAATGACTCAGGTTTCTGGCTCAAATCAAGTAT
CCAAGAGGCAACAGATGGACCCACTGCCTAGAATGAAGGAAAAGTCCGGTAAAACTGGAAGACACAATTC
CATGAATCAGAAGACAACAATGCCAAGAAAGGACCAGTTAAAAAGTTCCAAGATCAAACCATCTTCCA
AAACTCAACCCG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_177203.3](#), [NP_796177.2](#)

RefSeq Size: 2269 bp

RefSeq ORF: 645 bp


Locus ID: 320604

UniProt ID: [Q8BXX9](#)

Cytogenetics: 3 C

MW: 24.8 kDa

Product images:



Circular map for MR220920