

Product datasheet for **MR220985L3V**

Cfp (NM_008823) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cfp (NM_008823) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cfp
Synonyms:	BCFG; Pfc
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008823
ORF Size:	1395 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220985).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008823.3 , NP_032849.2
RefSeq Size:	1524 bp
RefSeq ORF:	1395 bp
Locus ID:	18636
UniProt ID:	P11680
Cytogenetics:	X 16.44 cM
Gene Summary:	A positive regulator of the alternate pathway of complement. It binds to and stabilizes the C3- and C5-convertase enzyme complexes.[UniProtKB/Swiss-Prot Function]



[View online »](#)