

Product datasheet for MR221000

Manf (NM_029103) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Manf (NM_029103) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Manf
Synonyms:	3230402M22Rik; AA407711; AA408789; AA673178; Armet; D18Mgi17; D17914
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR221000 representing NM_029103 Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTGGGCTACGCGGGCTGGCGGTAGCGCTGGCCCTGAGCGTGCTGCCTGACAGCCGGGCGCTGCGGC
CAGGAGACTGTGAAGTTTGTATTTCTTATCTGGGACGATTTTACCAGGACCTCAAAGACAGAGATGTCAC
ATTTTCACCAGCCACTATTGAAGAAGAAGTATAAAGTTTGGCGTGAAGCAAGAGGCAAAGAGAATCGG
TTGTGCTACTACATTGGAGCCACAGATGATGCTGCCACCAAGATCATCAATGAGGTGTCGAAGCCCCTGG
CCCACCATATCCCTGTGAAAAGATCTGTGAGAAGCTGAAGAAGAAAGACAGCCAGATCTGTGAAGTAAA
ATACGACAAGCAGATTGACCTGAGCACAGTGGACCTGAAGAAGCTCCGGGTGAAAGAGCTGAAGAAGATC
CTGGACGACTGGGGGAGATGTGCAAAGCTGTGCAGAAAAGTCTGACTATATCCGGAAGATAAATGAAC
TGATGCCTAAATACGCCCAAGGCAGCCAGCGCACGGACTGATCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR221000 representing NM_029103
 Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja1458_e02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_029103

ORF Size: 537 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_029103.4](#)

RefSeq Size: 2206 bp


RefSeq ORF: 540 bp

Locus ID: 74840

Cytogenetics: 9 57.98 cM

MW: 20.8 kDa

Product images:



Circular map for MR221000