

Product datasheet for **MR221047**

Nfasc (NM_001160316) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nfasc (NM_001160316) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Nfasc
Synonyms: AA387016; D430023G06Rik; mKIAA0756; NF
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR221047 representing NM_001160316
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCCAGGCAGCAGGCGCCACCCTGGGTCCACATAGCCCTCATCCTCTTCTCCTCAGCCTCGGAGGGG
CCATCGAGATCCGATGGACCTGACCCAACCCCAACTATCACCAAGCAGTCCGTGAAGGACCACATCGT
GGACCCTCGAGATAACATCCTGATTGAATGTGAAGCTAAAGGCAACCCCGCCCCAGTTTCACTGGACT
CGCAACAGCAGATTCTTCAACATTGCCAAGGACCCACGGGTGTCCATGAGGAGGAGATCTGGACCTTGG
TGATCGACTTCCGAGTGGTGGCGGCCTGAGGAATACGAAGGGGAGTACCAGTGCCTTGGCCCGAACA
ATTTGGCACTGCCTTAGCAACCGCATCCGCTGCAGGTGCCAAATCTCCCTGTGGCCCAAGGAAAAAC
CTAGACCCCGTCGTGGTTCAAGAGGGTGGCCCTTGACACTACAGTGCAACCCCCACCCGGCCTCCCGT
CCCCGTGATCTTCTGGATGAGCAGTCCATGGAGCCATCACCCAGGACAAGCGTGTCTCCAGGGTCA
CAACGGGACCTGTACTTCTCCAACGTCATGCTGCAGGACATGCAGACCGACTACAGTGCACCGCGC
TTTCACTTCAACCACACCATTGAGCAGAAGAACCCTTCAACCTCAAGTCTCACCAACAACCCCTATA
ATGACTCGTCTTAAGAAACCACCTGACATATATAGTGCCGAGGAGTTGCAGAAAGAAGCAGCCAGCTT
CATGTATCCCCAGGGCAGATCGAGCAGTCAAGTGGTTCTCCGTGGCATGGACCTGCTGCTGGAATGCAT
GCCTCTGGCGTCCAACACCAGACATTGCATGGTACAAGAAAGGTGGGGACCTCCCATCTAACAAGGCCA
AGTTCGAGAACTTTAATAAGGCTCTGCGCATCACCAATGTCTCTGAAGAGGACTCTGGGAGTATTTCTG
CCTGGCCTCCAACAAGATGGGCAGCATCCGGCACAGTCTCGGTGAGAGTAAAGGCTGCTCCATACTGG
CTGGATGAGCCCAAGAACCTGATCCTGGCTCCTGGGAAGATGGGAGGCTGGTATGCCGAGCCAATGGGA
ACCCGAAGCCGACCGTGCAGTGGATGGTGAATGGAGAGCCTTTACAATCGGCACCACCAATCCCAACCG
TGAGGTAGCTGGAGACACTATCATCTCCGGGATACTCAGATCAGCAGCAGGGCAGTGTACCAATGTAAT
ACATCCAATGAACATGGTACCTGCTGGCCAATGCCTTCGTGAGCGTGTAGATGTACCCCTCGGATGC
TGCTGCCCAGCAACAGCTCATCAGGGTATCCTTTATAACCCGGACACGGCTGGACTGTCCGTTCTTTGG
GTCTCCATCCCAACACTCCGATGGTTTAAAGATGGGCAAGGAAGCAACCTGGATGGCGGTAACCTACCAC
GTCTACGAAAACGGCAGTCTAGAAATCAAGATGATTGCAAGAGGACCAAGGCATCTACACCTGTGTGG



[View online >](#)

CCACCAACATCCTGGGCAAAGCCGAAAATCAAGTCCGCCTGGAGGTCAAAGACCCACCAGGATCTACAG
GATGCCCGAGGACCAGGTGGCCAAGAGGGGCACCACGGTGCAGCTGGAGTGCCGCGTGAAACATGACCCC
TCCTTGAAGCTCACAGTCTCCTGGCTGAAGGACGATGAGCCACTTACATTGAAACAGGATGAAGAAGG
AAGATGACTCCCTGACGATCTTCGGAGTGGCAGAGCGGGACCAGGGCAGTTACACGTGTATGGCCAGCAC
CGAGCTGGACCAGGACCTGGCAAAGGCCTACCTCACTGTTCTAGCTGATCAGGCCACTCCAATAACCGT
TTGGCTGCCCTACCCAAAGGGCGACCAGACCCAGGGACCTGGAGCTCACTGACCTGGCTGAGAGGA
GTGTGAGGCTGACCTGGATCCCAGGGGATGACAACAACAGCCCTATCACAGACTACGTGTTTCAGTTTGA
AGAGGACCAGTTCCAACCAGGGGTGGCATGACCACTCCAGGTTCCCAGGCAGCGTCAACTCAGCCGTC
CTCCATCTGTCCCATATGTCAACTACCAATTCAGAGTCACTGCTGTCAACGAGGTTGGGAGCAGCCACC
CCAGCCTTCCATCCGAGCGGTACCGAACAGTGGGGCACCCCTGAATCTAATCCCAGTGATGTGAAGGG
CGAAGGGACAAGAAAGAACAATATGGAGATCACGTGGACGCTATGAATGCTACCTTGCCTTTGGCCCC
AACCTACGCTACATTGTCAAGTGGCGACGGAGAGAGACCCGAGAGACTTGAACAATGTCACAGTGTGGG
GCTCTCGTACGTGGTGGGCAGACGCTGTCTACGTTCCCTATGAGATCCGAGTCCAGGCTGAAAATGA
CTTTGGGAAAGGCCCGAGCCTGACACCATCATTGGGTAATCCGAGAGATTATCCCAGGGCTGCGCCC
ACTGAAGTAAAAATCCGAGTCTGAACAGCACAGCCATCAGCCTTCAGTGAACCGAGTCTACTCTGACA
CGGTCCAGGGCCAGCTCAGAGAGTATCGAGCTTACTACTGGAGGGAAAGCAGTTTGTGAAAGAACCTGTG
GGTGTCTCAGAAGAGACAGCAGGCCAGCTTCCCTGGTGACCGTCCCCGGGGCGTGGTGGCCCGCCTGTT
CCCTACAGTAACTACAAGCTGGAGATGGTGGTGGTCAATGGGAGAGGTGACGGGCCCTCGAAGTGAACCA
AGGAATTCACCACCCAGAAGGAGTACCAGTGCCTCCAGGCGTTTCCAGAGTCCGACAGCCCAACCTGGA
GACCATCAACCTGGAGTGGGACCACCCAGAGCACCCCAACGGAATCCTGATTGGATACATCCTCAGATAC
GTGCCCTTTAATGGAACCAAACTGGGAAAGCAGATGGTGGAAAATTTCTCTCCAATCAGACCAAGTTCT
CTGTGCAGAGAGCAGACCAGTGTGCGGTTACCGCTTCTCCCTCAGTGCCAGGACACAGGTGGGCTCTGG
AGAAGCAGCCACAGAGGAGTCCCAGCACCTCCAATGAAGCTACTCCAAGTGCAGCTTACACCAATAAC
CAGGCAGACATCGCCACCCAGGGCTGGTTCATCGGGCTCATGTGTGCCATTGCCCTTCTGGTGTGATCC
TTCTCATCGTCTGCTTCAAGAGGAGTCGAGGTGGCAAGTACCAGTGCGGGAAAAGAAGGATGTCCC
CTTGGGCTCTGAAGACCCCAAAGAAGAAGATGGCTCATTTGACTACAGTGTGAGGACAACAAGCCCTG
CAGGGCAGCCAGACATCTCTGGATGGCACCATCAAGCAGCAGGAGAGCGATGACAGCCTGGTGGACTATG
GCGAAGGCGGCGAGGGCCAGTTCAATGAAGATGGCTCCTTTATTGGCCAGTACACTGTCAAAAAGGACAA
GGAGGAAACGGAGGGCAATGAGAGCTCAGAGGCCACATCACCAGTCAATGCCATCTATTCCTTGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >MR221047 representing NM_001160316
Red=Cloning site Green=Tags(s)


MARQQAPPWVHIALILFLLSLGGAIEIPMDLTQPPTITKQSVKDHI V D P R D N I L I E C E A K G N P A P S F H W T
 R N S R F F N I A K D P R V S M R R R S G T L V I D F R S G R P E E Y E G E Y Q C F A R N K F G T A L S N R I R L Q V S K S P L W P K E N
 L D P V V V Q E G A P L T L Q C N P P P G L P S P V I F W M S S S M E P I T Q D K R V S Q G H N G D L Y F S N V M L Q D M Q T D Y S C N A R
 F H F T H T I Q Q K N P F T L K V L T N N P Y N D S S L R N H P D I Y S A R G V A E R T P S F M Y P Q G T S S S Q M V L R G M D L L L E C I
 A S G V P T P D I A W Y K K G D L P S N K A K F E N F N K A L R I T N V S E E D S G E Y F C L A S N K M G S I R H T I S V R V K A A P Y W
 L D E P K N L I L A P G E D G R L V C R A N G N P K P T V Q W M V N G E P L Q S A P P N P N R E V A G D T I I F R D T Q I S S R A V Y Q C N
 T S N E H G Y L L A N A F V S V L D V P P R M L S A R N Q L I R V I L Y N R T R L D C P F F G S P I P T L R W F K N G Q G S N L D G G N Y H
 V Y E N G S L E I K M I R K E D Q G I Y T C V A T N I L G K A E N Q V R L E V K D P T R I Y R M P E D Q V A K R G T T V Q L E C R V K H D P
 S L K L T V S W L K D D E P L Y I G N R M K K E D D S L T I F G V A E R D Q G S Y T C M A S T E L D Q D L A K A Y L T V L A D Q A T P T N R
 L A A L P K G R P D R P R D L E L T D L A E R S V R L T W I P G D D N N S P I T D Y V V Q F E E D Q F Q P G V W H D S R F P G S V N S A V
 L H L S P Y V N Y Q F R V I A V N E V G S S H P S L P S E R Y R T S G A P P E S N P S D V K G E G T R K N N M E I T W T P M N A T S A F G P
 N L R Y I V K W R R R E T R E T W N N V T V W G S R Y V V G Q T P V Y V P Y E I R V Q A E N D F G K G P E P D T I I G Y S G E D Y P R A A P
 T E V K I R V L N S T A I S L Q W N R V S D T V Q G Q L R E Y R A Y W R E S S L L K N L W V S Q K R Q Q A S F P G D R P R G V A R L F
 P Y S N Y K L E M V V V N G R G D G P R S E T K E F T T P E G V P S A P R R F R V R Q P N L E T I N L E W D H P E H P N G I L I G Y I L R Y
 V P F N G T K L G K Q M V E N F S P N Q T K F S V Q R A D P V S R Y R F S L S A R T Q V G S G E A A T E E S P A P P N E A T P T A A Y T N N
 Q A D I A T Q G W F I G L M C A I A L L V L I L L I V C F I K R S R G G K Y P V R E K K D V P L G P E D P K E E D G S F D Y S D E D N K P L
 Q G S Q T S L D G T I K Q Q E S D D S L V D Y G E G G E G Q F N E D G S F I G Q Y T V K K D K E E T E G N E S S E A T S P V N A I Y S L A

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001160316

ORF Size: 3567 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001160316.1](#), [NP_001153788.1](#)

RefSeq Size: 9406 bp

RefSeq ORF: 3570 bp

Locus ID: 269116

Cytogenetics: 1 57.42 cM

MW: 134.2 kDa

Gene Summary: This gene encodes an L1 family immunoglobulin cell adhesion molecule with multiple IGcam and fibronectin domains. The protein functions in neurite outgrowth, neurite fasciculation, and organization of the axon initial segment (AIS) and nodes of Ranvier on axons during early development. Both the AIS and nodes of Ranvier contain high densities of voltage-gated Na⁺ (Nav) channels which are clustered by interactions with cytoskeletal and scaffolding proteins including this protein, gliomedin, ankyrin 3 (ankyrin-G), and betaIV spectrin. This protein links the AIS extracellular matrix to the intracellular cytoskeleton. This gene undergoes extensive alternative splicing, and the full-length nature of some variants has not been determined. [provided by RefSeq, May 2009]