

## Product datasheet for **MR221049L3V**

### Sema4f (NM\_001113481) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sema4f (NM_001113481) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sema4f
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001113481
ORF Size:	1662 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221049).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001113481.1</a> , <a href="#">NP_001106953.1</a>
RefSeq Size:	2735 bp
RefSeq ORF:	1665 bp
Locus ID:	20355
Cytogenetics:	6 35.94 cM



[View online »](#)

**Gene Summary:**

This gene encodes a member of semaphorin family of membrane-bound and secreted proteins that are involved in guiding axonal growth. The encoded protein is a transmembrane protein localized to the glutamatergic synapses via its association with a synapse-associated scaffolding protein. In oligodendrocyte precursor cells, the encoded protein contributes to the outward migration and differentiation. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, May 2015]