

Product datasheet for **MR221087L4V**

Asic2 (NM_001034013) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Asic2 (NM_001034013) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Asic2
Synonyms:	Accn1; ACIC2; BNaC1; BNaC1a; BNC1; Mdeg
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001034013
ORF Size:	1536 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221087).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001034013.2 , NP_001029185.1
RefSeq Size:	2643 bp
RefSeq ORF:	1539 bp
Locus ID:	11418
UniProt ID:	Q925H0
Cytogenetics:	11 48.43 cM



[View online »](#)

Gene Summary:

Cation channel with high affinity for sodium, which is gated by extracellular protons and inhibited by the diuretic amiloride. Also permeable for Li(+) and K(+). Generates a biphasic current with a fast inactivating and a slow sustained phase. Heteromeric channel assembly seems to modulate.[UniProtKB/Swiss-Prot Function]