

## Product datasheet for **MR221090L4V**

### Neu2 (NM\_001160163) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Neu2 (NM_001160163) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Neu2
Synonyms:	MBS; MSS; MTS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001160163
ORF Size:	1179 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221090).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001160163.1</a> , <a href="#">NP_001153635.1</a>
RefSeq Size:	1722 bp
RefSeq ORF:	1182 bp
Locus ID:	23956
UniProt ID:	<a href="#">Q9JMH3</a>
Cytogenetics:	1 D
Gene Summary:	Catalyzes the removal of sialic acid (N-acetylneuraminic acid) moities from glycoproteins, oligosaccharides and gangliosides.[UniProtKB/Swiss-Prot Function]



[View online »](#)