

Product datasheet for **MR221105**

Rnf157 (NM_027258) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rnf157 (NM_027258) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rnf157
Synonyms:	2610036E23Rik; A130073L17Rik; mKIAA1917
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>MR221105 representing NM_027258
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGGGGGCCCTGACGAGCCGGCAGCACGCGGGCGTGGAGGAGGTGGACATCCCGTCCAATTCCGTGTACC
 GCTACCCGCCCAAGTCCGGCAGCTATTTTGCCAGCCACTTCATCATGGGAGGAGAGAAGTTTGACTCCAC
 GCATCCGGAGGGTTACCTGTTTGCGGAGAACAGCGACCTGAACCTCCTGGGAAACAGGCCGGTGGCGTTC
 CCTTACGCTGCCCGCCTCCCAAGAACCTGTGAAGACTCTGAGAAGCCTGATCAACATCCGGAAGGACA
 CGCTGAGACTTGTCAAATGTGCTGAGGAAGTGAAGAGCCACGGAGAAGAGGCAGGCAAAGCTAAGGTCCA
 CTACAATGTGCGATTACCTTTGACACAGATGCCCGGTAGCCATCACCATCTACTACCAGGCCACTGAG
 GAATTCAGAATGGCATTGCCAGCTACATTCCAAAAGACAACAGCCTGCAGTCAGAGACCGTGCCTACA
 AGCGAGGGGTGTGCCAGAGTTTGCCTGCCCTCCCATACCGTGGACCCCTCGGAGTGGGCCGAAGAGGA
 GCTTGGCTTTGATCTGGACCGTGAAGTCTACCCTCTGGTAGTGCACGCTGTGGTGGATGAAGGAGATGAG
 TATTTTGGCCACTGCCACGTGCTTCTGGGCACTTTTGAGAAGCACCCGGACGGGACCTTCTGTGTGAAGC
 CCCTCAAGCAGAAGCAAGTGGTGGATGGGTGACGTACCTTCTCAGGAGATCTATGGGATTGAAAAACA
 GTACAACACACAAGATTCGAAGGTGGCAGAGGACGACGTGAGTGAACAACAGCGCTGAGTGTGTGGTGTGT
 CTCTCTGACGTCCGCGATACTTTGATCTGCCCTGCCGGCACCTCTGCCTCTGCAACACCTGCGCAGACA
 CCCTGCGCTACCAGGCCAACAAGTGCCTCATCTGCCGGTGCCTTCCGGGCACTGCTTCAGATCCGAGC
 CATGAGGAAAAAACTGGGCCCTCTGTCCCAAGCAGCTTAACCCCATCATCTCTCCAGACCTCGGAC
 TCTGAGGAGCATTATCTCAGAGAATCCCCCGGGCTATGAAGTGGTGTCTCTCCTCAGAGCCCTCA
 ACGGGCCCTCACCTCGTCCCCGGCAGTGCCTCCCTTACGTTCTTGGGGATGGCCACCTCTCAGGAAT
 GCTCCCGTCCTATGGCAGTGATGGCTACCTGCCCCCTGTTAGGACACTGTCCCCCTTGACCGCCTGTCT
 GACTGCAACAACCAAGGACTCAAACTCAAAAAGAGTCTATCCAAATCCATTTCCAGAAATCTTCTGTGC
 TTCACGAAGAGGAAGATGAGCGCTCTTGACGTGAGTCAAGACTCAGCTCTCTCAGAGGCTCTCAGTCCA
 GCACCCTGAAGAGGGACCCGACGTGACTCCAGAGAGTGAACCTCACCTGTCTCCTCAGGGGCTGTT
 GACCAGTCGTCTTGACAGGGACTCCGCTCTTCTACCATCTCTCCCCAGAAGACCCAGCCAGCAGCA
 GCCTGGCCAGTCAGTCATGTCCATGGCTTCTCCAGATCAGCACCGACACGGTGTCTCCATGTCTGG
 CTCCTACATCGCACCTGGTACAGAAGAAGAAGGAGAAGCCCTGCCTTCCCCCGAGCTGCTAGCAGGGCC
 CCTTCAGAAGGAGAGGAGACACCAGCAGAGTCCCAGACAGCAATTCGCTGGCCTTCCAGCTGGAGAGC
 AGGATGCAGAGGGAAATGATATCATAGAGGAAGAGGACAGATCACCTGTGCGAGAAGATGGCCAGAGGAC
 ATGCGCATTTCTAGGCATGGAGTGTGACAATAACAATGACTTTGACGTCGCAAGCGTGAAAGCACTGGAC
 AATAAGCTGTGCTCTGAGGTCTGCTTACCCGGTACCTGGCAGCATGAGGATAAACTGTCAACTGTGCGC
 ATACCCAGCGCCGGCGCTTATCTCCAGCAGCTGGAGGACCCTGAGGAGAACAGGCCTGTGTATGGGA
 TCCTATGGCTGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR221105 representing NM_027258
Red=Cloning site Green=Tags(s)


MGALTSRQHAGVEEVDIPNSNVYRYPKSGSYFASHFIMGGEKFDSTHPEGYLFGENSDLNFLGNRPVAF
PYAAPPPQEPVKTLRSLINIRKDTLRLVKCAEEVKSHGEEAGKAKVHYNVEFTFDTDARVAITYYQATE
EFQNGIASYIPKDNSLQSETVHYKRGVCQQFCLPSHTVDPSEWAEELGFDLDREVYPLVVHAVVDEGDE
YFGHCHVLLGTFEKHPDGTFCVKPLKQKQVVDGVSILLQEIYGIENKYNTQDSKVAEDDVSNSAECVVC
LSDVRDTLILPCRHLCLCNTCADTLRYQANNCIPCRLPFRALLQIRAMRKKLGPLSPSSFNPIISSQTS
SEEHSSSENIPPGEYVVSLLLEALNGPLTSSPAVPPLHVLGDGHLGMLPSYSGDGYLPPVRTLSPLDRLS
DCNNQGLKLLKSLSKSISQNSSVLHEEEDERSCSESDTQLSQRLSVQHPEEGPDVTPSENLTLSSGAV
DQSSCTGTPLSSTISSPEDPASSSLAQSVMSMASSQISTDTVSSMSGSYIAPGTEEEGEALPSPRAASRA
PSEGEETPAESPDSNFAGLPAGEQDAEGNDIIEEDRSPVREDGQRTCAFLGMECDNNNDFDVASVKALD
NKLCSVEVCLPGTWOHEDNTVNCRHTORRRLLSSSLEDPEENRPCVWDPMAY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja2207_h05.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_027258

ORF Size: 2043 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_027258.1](#), [NP_081534.1](#)

RefSeq Size: 4655 bp


RefSeq ORF: 2046 bp

Locus ID: 217340

Cytogenetics: 11 E2

MW: 74.9 kDa

Product images:



Circular map for MR221105