

Product datasheet for **MR221107**

Mum111 (NM_175541) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mum111 (NM_175541) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Mum111
Synonyms:	9430086K21Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR221107 representing NM_175541
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGATGGCAAATATGTCCTGTGTAATTGGAAAGACCAATTATGGCCAGCAAAGTTTTAGATAGGTCTG
 AAAGTCCATCAGAGAGTAAGAGGAAAAATACATCGTCCCTAGAAGTTGAAATACTTTCACTGGATGAGAA
 AATTACCGTGAAAGCACAGACACAAAAGTCTAAGCAAATCTGCTGTTGAAGCCATTATGTCCTCTCTA
 GCAGTACAGTCAGAGGTCATCATTGCACCTAGAGAGGAGACAGCCTATGAAAGATCACTGAAAATGGCAC
 TGGAAATGGTGAAGAAGGAACCAATCTCAGCCAGGAAAGCATGTCAGAAGAACAACCTACAGCTACAGC
 GTCTGAAAATGTGCCAGAACAGCCACCTGATTCACCTCCTCATAAAAAGTTCCGCAAACCTGAAAGCAAC
 ACCCAGGAGGACTCTGCTTCCATACTGCTGTGCTCAGAGAGTGATGATTCCATGACCGATGATAAGTTGC
 AGGTGCACACAACCGGTGAGAGCATGCCAAGTGAGATGGACACAAAGGCAACAGAAAACCTTAGGCTGTTG
 CCAAACAGACCCTTCTTGGCAGATGAGGATGATAAAAAGGAAGAAAAGAAAAGATTGACATCTCGGCA
 ATCATGTCTGTGAATTTATCACTCAAAGAAGAAAGCGAATATATTAAGAAGAGAAAATTCGTCCCTCAT
 CAGAAGATCTCACTGTCCCAAAGAGGAGTCACAAGACATCCTCCAGAAGCGCCTCTGGCTGTTTCCTC
 TGAATGTTCTATTGTCTCAGAGAATAACATGGAAGATCCTGGTGAGGGCCCATCAAATCAGAATCTAGTC
 TCCTATGCCAACCAAAATCAGTCTTCTGTGGAATCCGATGTAGGTGCCGAGACATCCACTGCAGGGTGT
 CAGGGGACTTTCAGGTTTCCCTTCTACCCGTGATACAGTCAGCAGTGATCTACTGCTCCAGAGACTGGA
 TCTAGAAGATCTTGAGGAAGAAGCCGAGCTTCAGGCAAGCTTTTGTCTTTAAATCCTGCCAGTGCCGCT
 GCGTTAGAAAATGATAATGAAGATGATGATGAAGACCTTCCACGTTTTCATCTCCATTATGAGACACGTG
 CATTGAAAATGGAATGATAGTGTGGTTAAATATCAAAAATACCCATTTTGCCAGCAGTGATCAAAAAG
 CATTAGGCGGAAAAGAGAGGAAAAGCGAGTGTGCTTTTGGTTGAGGCAGACATGAGTCCTCAAAGAAAGGC
 GTTAGAGTATCTTTGAGGAGGCTGAAGAAATACGACTGTAAAGAGAAAACAGGCACTAGTGAGAAAAGCCA
 GGGAGGAGTACCGGGAGAGTATCGATTGGTGCCTGCTGACTGATTTGTGACTACCGAGTTAGACTAGGTTG
 TGGCTCTTTACCGGCTCGTCTTTGAGTATTATGCTGCTGACATCAGTTATCCAGTCAGGAAAATCATC
 AAACAAGATACCTTCAGAAAATATTTCCAAAGCTATACAATGAAAACGTGGGGGAGCAGTTGCCTATGG
 CATCCCAAGCCAAGAGAGTATCTTTCCAGAAAATCTCCCTGATCGGATGAAGCCTGCTCGGGACCGAGC
 TAACAAGAACCTGGTAGATTTTATTGTCAATGCAAAGGAACAGAGGACCACCTCTTGGCATTTTAAAG
 GGCACCAAAAAATCCAAGTGGCTGAAATCATTTCTCAATGCAAAGAGTTTACACCTGTATTGAAACCT
 ACTTTGAAGATGAAGATCAACTGGATGAGGTAGTGAAATACCTACAAGAAATTTATAAGCAAATTGACCA
 GAAGATGCTGACTCTGATAAAAGATGACAAAATTAAGTTTGTCTTGAAGTTCTTCTACCAGAAGCAATT
 ATTTGCTCAATTTCTGCTGTTGATGGCTTAGATTATGAGGCAGCTGAGGCAAAGTATCTAAAGGGACCAT
 CCCTCGGCTGTAGGGAGAGAGAATTATATGATTCCAAAATCCTATTTGAAAAGAGAAGGAGGTCATTACC
 AAACGAAGGCCAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR221107 representing NM_175541
Red=Cloning site Green=Tags(s)

```
MDGKYVLCNWKDQLWPAKVLDRSESPSESKRKNTSSLEVEILSLDEKITVESTDTKVLKSAVEAIMSSL
AVQSEVIIAPREETAYERSLKMALVMKEGTNL SQESMSEEQPTATASENVPEQPPDSPPHKKFRKLESN
TQEDSASILLCESDDSMTDDKLQVHTTGESMPSEMDTKATENLGCCQTDPSLADEDDKKEEKKIDISA
IMSVNLSLKEESEYIKEEKFPSSDEL TVPKEESQDILPEAPLAVSSECSIVSENMEDPGEGPSNQNLV
SYANQNQSSVESDVGAETSTAGCSGDFQVSLPTRDTVSSDLLLQRLDLEDLEEEARASGKLLSLNPASAA
ALENDNEDEDLPRFILHYETRAFETGMIVWFKYQKYPFWPAVTKSIRRKERKASVLLVEADMSPQKKG
VRVSLRRLKKYDCKEKQALVEKAREEYRESIDWCVSLICDYRVRLGCGSFTGSFFEYAAADISYPVRKII
KQDTRFNIFPKLYNENVGELPMASQAKRVSFQKILPDRMKPARDRANKNLVDFIVNAKGTEDHLLGILK
GTTKSKWLKSFNAKSFPCIEITYFEDEDQLDEVVKYLQEIKYQIDQKMLTLIKDDKIKFVLEVLPEAI
ICSI SAVDGLDYAAEAAYLKGPSLGCRRERELYDSKILFEKRRRSLPNEG
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_175541

ORF Size: 2043 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_175541.3](#), [NP_780750.2](#)

RefSeq Size: 4523 bp

RefSeq ORF: 2046 bp


Locus ID: 245631

UniProt ID: [Q4VA55](#)

Cytogenetics: X F1

MW: 77 kDa

Product images:



Circular map for MR221107