

## Product datasheet for **MR221112**

### **Mum111 (NM\_001164630) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Mum111 (NM_001164630) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Mum111
Synonyms:	9430086K21Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MR221112 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCCGCATCGCC

ATGGATGGCAAATATGTCCTGTGTAATTGGAAAGACCAATTATGGCCAGCAAAGTTTTAGATAGGTCTG  
AAAGTCCATCAGAGAGTAAGAGGAAAAATACATCGTCCCTAGAAGTTGAAATACTTTCACTGGATGAGAA  
AATTACCGTGGAAAGCACAGACACAAAAGTCTAAGCAAATCTGCTGTTGAAGCCATTATGTCCTCTCTA  
GCAGTACAGTCAGAGGTCATCATTGCACCTAGAGAGGAGACAGCCTATGAAAGATCACTGAAAATGGCAC  
TGGAAATGGTGAAGAAGGAACCAATCTCAGCCAGGAAAGCATGTGAGAAGAACAACCTACAGCTACAGC  
GTCTGAAAATGTCCAGAACAGCCACCTGATTCACCTCCTCATAAAAAGTTCGCAAACTTGAAAGCAAC  
ACCCAGGAGGACTCTGCTTCCATACTGCTGTGCTCAGAGAGTGATGATTCATGACCGATGATAAGTTGC  
AGGTGCACACAACCGGTGAGAGCATGCCAAGTGAGATGGACACAAAGGCAACAGAAAACCTTAGGCTGTTG  
CCAAACAGACCCTTCTTGGCAGATGAGGATGATAAAAAGGAAGAAAAGAAAAGATTGACATCTCGGCA  
ATCATGTCTGTGAATTTATCACTCAAAGAAGAAAGCGAATATATTAAGAAGAGAAAATTCGTCCCTCAT  
CAGAAGATCTCACTGTCCCAAAGAGGAGTCACAAGACATCTCCAGAAGCGCCTCTGGCTGTTTCCTC  
TGAATGTTCTATTGTCTCAGAGAATAACATGGAAGATCCTGGTGAGGGCCCATCAAATCAGAATCTAGTC  
TCCTATGCCAACCAAAATCAGTCTTCTGTGGAATCCGATGTAGGTGCCGAGACATCCACTGCAGGGTGT  
CAGGGGACTTTCAGGTTTCCCTTCTACCCGTGATACAGTCAGCAGTGATCTACTGCTCCAGAGACTGGA  
TCTAGAAGATCTTGAGGAAGAAGCCGAGCTTCAGGCAAGCTTTTGTCTTTAAATCCTGCCAGTGCCGCT  
GCGTTAGAAAATGATAATGAAGATGATGATGAAGACCTTCCACGTTTCATCTCCATTATGAGACACGTG  
CATTTGAAACTGGAATGATAGTGTGGTTAAATATCAAAAATACCCATTTTGCCAGCAGTGATCAAAAAG  
CATTAGGCGGAAAAGAGAGGAAAAGCGAGTGTGCTTTTGGTTGAGGCAGACATGAGTCCTCAAAGAAAGGC  
GTTAGAGTATCTTTGAGGAGGCTGAAGAAATACGACTGTAAAGAGAAAACAGGCACTAGTGAGAAAAGCCA  
GGGAGGAGTACCGGAGAGTATCGATTGGTGCCTGCTGACTGATTTGTGACTACCGAGTTAGACTAGGTTG  
TGGCTCTTTACCGGCTCGTCTTTGAGTATTATGCTGCTGACATCAGTTATCCAGTCAGGAAAATCATC  
AAACAAGTACCTTCAGAAAATATTTCCAAAGCTATACAATGAAAACGTGGGGAGCAGTTGCCTATGG  
CATCCCAAGCCAAGAGAGTATCTTTCCAGAAAATCTCCCTGATCGGATGAAGCCTGCTCGGGACCGAGC  
TAACAAGAACCTGGTAGATTTTATTGTCAATGCAAAGGAACAGAGGACCACCTCTTGGCATTTTAAAG  
GGCACCAAAAAATCCAAGTGGCTGAAATCATTTCTCAATGCAAAGAGTTTCACACCCTGTATTGAAACCT  
ACTTTGAAGATGAAGATCAACTGGATGAGGTAGTGAATACCTACAAGAAATTTATAAGCAAATTGACCA  
GAAGATGCTGACTCTGATAAAAGATGACAAAATTAAGTTTGTCTTGAAGTTCTTCTACCAGAAGCAATT  
ATTTGCTCAATTTCTGCTGTTGATGGCTTAGATTATGAGGCAGCTGAGGCAAAGTATCTAAAGGGACCAT  
CCCTCGGCTGTAGGGAGAGAGAATTATATGATTCAAAATCCTATTTGAAAAGAGAAGGAGGTCATTACC  
AAACGAAGGCCAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA


**Protein Sequence:** >MR221112 protein sequence  
Red=Cloning site Green=Tags(s)

```
MDGKYVLCNWKDQLWPAKVLDRSESPSESKRKNTSSLEVEILSLDEKITVESTDTKVLKSAVEAIMSSL
AVQSEVIIAPREETAYERSLKMALVMKEGTNL SQESMSEEQPTATASENVPEQPPDSPPHKKFRKLESN
TQEDSASILLCESDDSMTDDKLQVHTTGESMPSEMDTKATENLGCCQTDPSLADEDDKKEEKKIDISA
IMSVNLSLKEESEYIKEEKFPSSDEL TVPKEESQDILPEAPLAVSSECSIVSENMMEDPGEGPSNQNLV
SYANQNQSSVESDVGAETSTAGCSGDFQVSLPTRDTVSSDLLLQRLDLEDLEEEARASGKLLSLNPASAA
ALENDNEDDDEDLPRFILHYETRAFETGMIVWFKYQKYPFWPAVIKSIRRKERKASVLLVEADMSPQKKG
VRVSLRRLKKYDCKEKQALVEKAREEYRESIDWCVSLICDYRVRLGCGSFTGSFFEYAAADISYPVRKII
KQDTRFNIFPKLYNENVGELPMASQAKRVSFQKILPDRMKPARDRANKNLVDFIVNAKGTEDHLLGILK
GTTKSKWLKSFNAKSFPCIEITYFEDEDQLDEVVKYLQEIKYQIDQKMLTLIKDDKIKFVLEVLPEAI
ICSI SAVDGLDYAAEAAYLKGPSLGCRRERELYDSKILFEKRRRSLPNEG
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001164630

**ORF Size:** 2043 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001164630.1](#), [NP\\_001158102.1](#)

**RefSeq Size:** 4600 bp

**RefSeq ORF:** 2046 bp


**Locus ID:** 245631

**UniProt ID:** [Q4VA55](#)

**Cytogenetics:** X F1

**MW:** 77 kDa

**Product images:**



Circular map for MR221112