

Product datasheet for **MR221124L3V**

Blzf1 (NM_001160208) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Blzf1 (NM_001160208) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Blzf1
Synonyms:	1700030G05Rik; Blzf1l; Golgin-45; Jem-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001160208
ORF Size:	1212 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221124).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001160208.1 , NP_001153680.1
RefSeq Size:	3927 bp
RefSeq ORF:	1212 bp
Locus ID:	66352
UniProt ID:	Q8R2X8
Cytogenetics:	1 H2.2
Gene Summary:	Required for normal Golgi structure and for protein transport from the endoplasmic reticulum (ER) through the Golgi apparatus to the cell surface.[UniProtKB/Swiss-Prot Function]



[View online »](#)