

Product datasheet for **MR221177**

Zfp879 (NM_173387) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zfp879 (NM_173387) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Zfp879
Synonyms:	Znf879
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR221177 representing NM_173387
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGACCACGAGGTTGCTGCCAGCCAGGTTTCAGGAGCCTGTGACGTTTCAGGGATGTGGCCGTGTCCTTCA
 GCCAGGATGAGTGGTGGACCTGGACCTGCTCAGAGAACCTTGACCGAGAGGTGATGCTGGAGAACTA
 CAGCAACCTGGTCTCGCTGGGGATCCTGTTTTCAAGCCAAAGGTCATCGCTCAGTTAGAACAAGCAGAA
 GACTTCCCGATGGTGAAAGTGGGATGTTTCAAGGTGTGTACCTCGGATGGAAGAGCTTGTTTGAGCCCA
 TAGTTTTCAAAGAAGAAAATAAGAAAGTCTTGAACCAAGAGAGTGGATAAATTTTTTAACTTCATATT
 GGGAAAAACGTGTGTAATGAGGAACAGTTAGAGGAGCAACAGGCAATAAGAAAAACCTTTTCAGCGAG
 ATGCTGCTTTCTGTCAGAAGAGGCTGTAGACAAGAGCGAGGCTTCCAGGGGTGGACTCGTGAGAAATC
 CCGGCATAAAGCTCCCATTTCTAGGAAACCTGGAGTTGCCGCCGAGGGGAAAGCCCCACTTACAGCA
 GTACTCAGTTCTGTTCAAACAACCTGGGATTCAACCGAGCCGCCAGCTTTGTAATGCAACATCTGTGGG
 AAAGTCTTCTTCCACAGTTCTCTCCCTGAGCAAGCATCAGAGGACTCACACTGGAGAGAAGCTCTACAAAT
 GTCAGGGCTGTGCGAAAGCTTTAGCCAACGCTCATCCCTGGCACAGCACCTGAGAGTCCATACTGGGGA
 AAAGCCGTACCTCTGTAGTGACTGTGGGAAGGCCTTCAGCTTCACCACCTCTCTCATTGGGCATCAGAGG
 ATGCACACTGGAGAGAGACCTATGAATGTAAGGAATGTGGCAAGACCTTCAAAGGCAGCTCCTCCCTTC
 ATAATCACCAACGGATCCACACTGGAGAAAAGCCTTATAAGTGTAAAGAGTGTGGGAGAGCCTTTAGCCA
 GTGCTCATCTCTGATTGAGCACCATAGGATTACACTGGGGAGAAACCATATGAATGCAGCCAGTGTGGG
 AAGGCCTTTACATCAATATCACGGCTGAGCAGACACCACCGCTCCATACCGGAGAGAAGCCCTTCACT
 GTAATGTGTGTGGCAAAGTGTTCAGTTACCCTCGGCTTGACCATCCATCAGAGAACTCACACTGGAGA
 GAAACCATATGCATGCAAGGAATGTGGGAAGCCTTCAGCCAAAGCTCAGCACTTACCCAACATCAAAGA
 ATCCACACTGGAGAAAAGCCGTACAAGTGTGCCGAATGCGGAAAGGCCTTTTCTGGCTCTCACGGCTTA
 ACATCCACCACAGGATCCACACGGGGGAGAAACCGTATCACTGTAAAGAGTGTGGGAAAGCCTTCAGTTC
 CCACTCTGCGGTCAATACGCACAGAAAGATCCACACTGGTGAGAAACCTTATAATGCAGTGACTGTGAA
 AAAGCCTTCAATCAAAGCTCGGCCCTTATTCAGCATCAGAGAATCCACACCGGAGAGAAGCCGTTTAACT
 GTAAAGTTTGTGGGAAAGCCTTCGACAGAGTTCTCTCGCTTATGACACATATGAGGATCCACACAGGAGA
 AAGGCCGTACAGATGCGAAGCTTGGGGAAAGCGTTCAGTCAGAGCTCATCTCTTGCCAACCACAGAAG
 ACACATTAC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR221177 representing NM_173387
 Red=Cloning site Green=Tags(s)


MTTRLLPAQVQEPVTRFDVAVSFQDEWLHLDPAQRTL YREVMLENYSNLVSLGILFSKPKVIAQLEQAE
 DFPMVESGMFQGVYLGWKSLEFPIVSKEENKKVLKNQRVDKFFNFILGKTCVNEEQLEEQPGNKNPFSE
 MLLSVRRGCRQERGFPGVGLVRNPGIKLPFRKPGVAAGGGKPHLQQYSVLFKQLGFNRARQLCKCNICG
 KVFFHSSSLSKHQRTHTGEKLYKQGCRAFSQRSSLAQHLLRVHTGEKPYLCSDCGKAFSFTTSLIGHQR
 MHTGERPYECKEKGTFKGSSSLHNHQRHTGEKPYKCNECGRAFSQCSSLIQHHRHTGEKPYECSQCG
 KAFTSISRLSRHHRVHTGEKPFHCNVCGKVF SYHSALT I HQRTHTGEKPYACKECKGKAFSQSSALTQHQR
 IHTGEKPYKCAECGKAFSWL SRLNIHHRHTGEKPYHCKECKGAFSSHSVNTHRKIHTGEKPYKSDCE
 KAFNQSSALI QHQRIHTGEKPFNCKVCGKA FRQSSSLMTHMRIHTGERPYRCEACGKAFSQSSSLANHQK
 THY

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI


Cloning Scheme:



ACCN:	NM_173387
ORF Size:	1689 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_173387.3 , NP_775563.1
RefSeq Size:	2427 bp
RefSeq ORF:	1692 bp
Locus ID:	214779

UniProt ID: Q8BI99
 Cytogenetics: 11 B1.3
 MW: 64.1 kDa
 Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR221177