

Product datasheet for **MR221182**

Zscan22 (NM_001001447) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zscan22 (NM_001001447) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Zscan22
Synonyms:	AW227621; D530006B18Rik; Hkr2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR221182 representing NM_001001447
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCCATCCCCAAGAGCCCTCTGAGCCAGTACTCTGGGAACAGGATAGCTTCCTTCAGGTGAAGGTGA
 AAGAGGAGGAAGAAGCTAGCGACTCCCCGAGCCAGGAATCCAGCCCTAGCTTTACTGCCACCCTGAGGC
 TGCACGATTGCGCTTCAGGCACTTCGCTATGAGGAGGCATCCAGCCCTCATCAGGCGCTAGCCAACTG
 CGAGAGCTATGTTGCCAGTGGCTGAGGCCAGAGTCGTCCTCAAAGAACAGATGCTGGAGTTGCTGGTGC
 TGGAGCAGTTCCTGGGTGCGCTGCCCTGAGATCCAGGCCTGGTGGGTGCTCAGTGCCCAAAGAGTGG
 AAAGGAGGCTGCCGTTCTGGTGGAGGATATGACTCAACTGCTGGACAGGAGAGGATGGGAGCCAGGAGTG
 GAACGTGAAGAGGCAAGTTGCAAACAGAGCAATACAGATGAATTGGAGCCACCAAAAATGGCGACTGAAA
 CTGTCATGGGAAGCGTTTGC AAAATCTACCCTTGCCACACTTGTAACCTGAGAGCCACTCAGAGAG
 CCAGCCAGAATTCCTAGGAGCACTCTGGATGAAGTCTACTGCCCAAGAGATGGATTTAGAAAAGCTCTG
 GGGCCTCATATGGATGCCCCAAAGACCAGCCTGGCCATGAGTCTAACACCTCAGGAAACGGTTCCAACA
 TGTGGCCCAACTCCCTCCCAAGACAAGGCTTCTCAGAGGAGAAATTTGGCCATTACTTGATAATGA
 GACAGTGCCTCCAGATACTTGCTCAGAGAAGAAATCCTCCAAGGATAGTGAATGTTTGAAAACATCCAG
 AACACTTCAGCTCTGGAGGCCACCAGAAGAGCCACTCTCAGAAGACGCCCTATGCCTGTACTGAGTGTG
 GGAAGGTATTCAGTCGGAGCACCCACCTTGCCAAACACAGGTTGTTACACAGGAGCAAAGCCCATGC
 GTGTAAGAGTGTGGTAAAGCCTCAGTCGCGTTGCCACCTGACACAGCACCTGAGGATTCATACTGGG
 GAGAAGCCTTATAATGTGAGGAATGTGGTAAAACCTTCAGCCGACGACCCACCTCAGCAGCACCAGC
 GAGTACACACAGGCGAGCGGCCCTACGAGTGTGATACTTGGGCAAAGCCTTCAGCCAGAGCACCACCT
 GACACAGCACAGCGCATCCACACCGGGGAGAAGCCCTACAGGTGTGATGTGTGTGGAAAAGCCTTCAGT
 GACTGTTGAGCTCTGGTCCGGCACCTGAGAGTCCACTCTGGAGAGAAGCCTTATCAATGTAAAGACTGCC
 CGAAGGCCTTCGCACAAAGCTCCTCCCTCATTGAGCACCAGCGCACTCACACAGGTGAAAAGCCTTATAA
 ATGCAAGTACTGTGGAAAGCTTTAGCCGAGCTCAGCTCTCATGGTGCACCTAAAGATCCATATCACC
 GTGACTCGAAGCACCCCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR221182 representing NM_001001447
 Red=Cloning site Green=Tags(s)

MAIPKSPLSPVLWEQDSFLQVKVEEEEASDSPSQESSPSFTHPEAARLRFHFYEEASSPHQALQ
 RELCCQWLRPSSSKEQMLELLVLEQFLGALPPEIQAWVGAQCPKSGKEAAVLVEDMTQLLDRRGWEPGV
 EREEASCKQSNTDELEPPKMATEVMGSVLPKSTLAHTCKPESHSESQPEFLGALWMKSTAQEMDFRKAL
 GPHMDAPKDQPGHESNTSGNGSNMWPFPSPQDKASSEKFGPLLDNETVPPDTCSEKSSKDSECLKTFQ
 NTSALEAHQKSHSQKTPYACTECGKVF SRSTHLVQHVVHTGAKPHACKECGKAFSRVAHLTQHLRIHTG
 EKPYKCEECGKTF SRSTHLTQHQRVHTGERPYECDTCGKAFSQSTHLTQHRIHTGEKPYRCDVCGKAFS
 DCSALVRHLRVHSGEKPYQCKDCPKAFAQSSSLIEHQRTHTGEKPYKSDCGKAFSRSSALMVHLKIHT
 VTRSTP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_001001447

ORF Size: 1488 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001001447.3](#), [NP_001001447.1](#)

RefSeq Size: 2620 bp


RefSeq ORF: 1491 bp

Locus ID: 232878

Cytogenetics: 7 A1

MW: 55.9 kDa

Product images:



Circular map for MR221182