


Product datasheet for MR221232L4

Itga10 (NM_001081053) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Itga10 (NM_001081053) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Itga10
Synonyms:	A630048L14
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221232).
Restriction Sites:	SgfI-RsrII
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_001081053
ORF Size:	3501 bp



[View online »](#)

