

## Product datasheet for **MR221253L4V**

### **Kcna10 (NM\_001081140) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Kcna10 (NM_001081140) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Kcna10
Synonyms:	Gm1962; Kcna8; Kv1.8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001081140
ORF Size:	1533 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221253).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001081140.1</a> , <a href="#">NP_001074609.1</a>
RefSeq Size:	1826 bp
RefSeq ORF:	1536 bp
Locus ID:	242151
UniProt ID:	<a href="#">B2RQA1</a>
Cytogenetics:	3 46.64 cM



[View online »](#)

**Gene Summary:**

Mediates voltage-dependent potassium ion permeability of excitable membranes. Assuming opened or closed conformations in response to the voltage difference across the membrane, the protein forms a potassium-selective channel through which potassium ions may pass in accordance with their electrochemical gradient. The channel activity is up-regulated by cAMP. [UniProtKB/Swiss-Prot Function]