

Product datasheet for MR221254

Adamts15 (NM_025629) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Adamts15 (NM_025629) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Adamts15
Synonyms: 2010109H09Rik; Thsd6
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR221254 representing NM_025629
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAATCACCTGACCCTGCTCTTCACTGCCTCTTGGTAACCCTGTGCTCCCTCAAACCCCGCCTCT
TCCCGCTTCTCCGTTCTTGCTGTGGACTCTCCTGAGCTGTGATTTGAGGGGTGCTGCCAGGGTCCGGG
AGAGTGGACGCCCTGGGGTCTTGGAGCAGGTGTCCAGCTCTTGTGGTCGGGGCTCTCAGTGCCGAGA
CGCCAATGTGTCCGACTTCCCGGAGAGGAACCTTGTGGGGCGACAGCCATGAATACCGGTGTGCCAGC
TGCCGGACTGCCCTCCAGGGACCACACCATTCGAGACCTGCAGTGTTCTCTACAATGCCACCCCTGT
CCTGGACACCCAGAAAACCTACCAGTGGGTGCCCTTCCACGGGGCGCCCAATGTGTGTGACCTCAACTGC
CTAGCTGAAGGCCACGCCTTCTACCACAGCTTTGGCCGCGTGTGGATGGCACCCCTGCACCCCGGTA
CCCAGGGTCTCTGCGTGGCAGGCCGCTGTCTGAGCGCAGGCTGTGACGGGATACTGGGCTCCGGCACCT
CGAGGACCGCTGCGGCCCTCTGTGGTGGCGCTAACGACTCTTGCCTGTTCTGTCAGCGCGTGTCCGCGAC
GCCGGTACTGTCCATTCTGACTCACCACAGCGTTGAGAAACACAGTGGAGAGGGGTGCCCTTTAA
CAAGTTGGCTCCTGGGAGCCAAGGCTCTTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR221254 representing NM_025629
Red=Cloning site Green=Tags(s)

MAIHPDPALHCLLVTLCSLKPRLLPFLPFWLLSCDLRGAQGPGEWTPWGSWSRCSSSSCGRGVSVRR
 RQCVRLPGEELCWGDSHEYRVCQLPDCPPGTTFFRDLQCSLYNGHPVLDTKQTYQWVFFHGAPNVCDLNC
 LAEGHAFYHSFGRVLDGTPCTPGTQGLCVAGRCLSAGCDGILGSGTLEDRCGLCGGANDSCLFVQVRFRD
 AGDCPILTHHQRLRNTVERGVPLTSLAPGSQGSC


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_025629
ORF Size: 732 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_025629.2 , NP_079905.1
RefSeq Size:	989 bp
RefSeq ORF:	735 bp
Locus ID:	66548
Cytogenetics:	10 C1
MW:	36.2 kDa