

## Product datasheet for MR221280

### Rhox2f (NM\_001085356) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rhox2f (NM\_001085356) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Rhox2f  
**Synonyms:** EG434764; Rhox2.6; Rhox2e  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR221280 representing NM\_001085356  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGGATCGCC

ATGGAGCGTCAAAGCATCAATTACCTGCTTGTGTGGACCTGAGGAGGATGAGGAAAAATGCGAATGGTG  
 TAAAGACTCTGATGGTCTTGCTGGCTGGAGACGGAAGAAACGAGGGAGAGAGTGGACGGGCTGCCTGG  
 GTCAGGAGCCTCAGCAGCGGAAGGATACAGAGCAGGAGAAATAAGTGCAGGTGGGCTGCTGCGCAAGTA  
 GCCGACCTCATGGATGACAGCAACCAAGAGGACCTTGGTGCCACTGCCTGTGACCAGGAGAAGGAGAAGC  
 AGCCAGAGGAGCCAGTCCCTGATTCCATGGGAGATTTGGAAAATGTAAGCGTGTGTCCGGGCCGTGGTC  
 CACTGTTAATCCTGTGAGAGTGTGGTGCCAAATCCGGCACCTTTGGCGACACAGCTCAATGTGCTG  
 CAACTACAAGAGCTGGAGAGCATCTCCAGTGAATCACTACATCAGCACTAAGGAGGCAAAATCGCCTGG  
 CAAGATCCATGGGTGTGAGTGAAGCCACAGTGCAGGAATGGTTTTTGAAGAGGAGAGAAAAATACAGGAG  
 TTATAAGAGGCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR221280 representing NM\_001085356  
 Red=Cloning site Green=Tags(s)

MERQSYNYLLDVGPEDDENANGVKTLMLVLLAGDGRNEGESGRGLPGSGASAAEGYRAGEISAGGPAAQV  
 ADLMDDSNQEDLGATACDQEKEKQPEEPVPDMSGDLENVSRVSGPSTVNPVVRVLPKFRHLWRHSFNVL  
 QLQELSIQCNHYISTKEANRLARSMGVSEATVQEWFLKRREKYRSYKRL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV




[View online »](#)

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9072\\_d12.zip](https://cdn.origene.com/chromatograms/mm9072_d12.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001085356

**ORF Size:** 573 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001085356.1](#), [NP\\_001078825.1](#)

**RefSeq Size:** 779 bp


**RefSeq ORF:** 576 bp

**Locus ID:** 434764

Cytogenetics: X A3.3

MW: 21.2 kDa

**Product images:**



Circular map for MR221280