

## Product datasheet for **MR221316L4V**

### Sult2a1 (NM\_00111296) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sult2a1 (NM_00111296) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sult2a1
Synonyms:	mSTa1; ST2A1; Std; Sth1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_00111296
ORF Size:	855 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221316).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_00111296.2</a> , <a href="#">NP_001104766.1</a>
RefSeq Size:	1021 bp
RefSeq ORF:	858 bp
Locus ID:	20859
UniProt ID:	<a href="#">P52843</a>
Cytogenetics:	7 8.02 cM
Gene Summary:	This is one of seven sulfotransferase family 2A genes in a chromosome 7 A1 cluster. [provided by RefSeq, May 2010]



[View online »](#)