

Product datasheet for **MR221330**

Lonrf2 (NM_001029878) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Lonrf2
Synonyms:	2900060P06Rik; AI851349
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence: >MR221330 representing NM_001029878
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAAGAACTATGATCAAGCCCTGCAGGATGCTGAAGCAGTTTGTCTAGAGGGAACCTCTCCTGACCAAGG
 GGCACCATATAAAAGCCAGATTCTTTCTGGATTGGGAAGACATAGAGAAGTGTAAAGGAATTTATCTA
 CTGCCTTGCTTTAAACCTGAGTCAACTCGGCCGAAGAAAGAAACAGAGAAGGTAATATGTGAGGTGTTT
 TTTTCAACTTCAGAAAGTGAGCACCACCACTCAACATCTTCCACTGAACTGGACCAGAGGCCCTGTGTG
 AAGGCCAGACGAATCCCAAGTATCTCTGGAAGAAGCAGGCGGTACGCTAATGCAGACAATCCTAAGAC
 TCCATCTGAGAAAGTCAGATGCACCCGCGGATATCAACTCTTCGGTTTTATATTTTATTCTGGGCCTCCAC
 TGTGAAGAGGATAAAAAGGCTCTGGAAGGCATCGTCCAGCGGCACCCAGCAGCACCTTGAAGAGGCAGC
 TTCCCAAGTGATGCCAGGACGATGAGGAGCTCAAAGCAATACCCCTGAGAAGATCCCAAGAAAGATGC
 AGACTCCCCACCTCAGAGGAATGCGAGCTCGTTAGAGGAAGAGCCAGAGTTACAATCGATGCGACTGAC
 TTTGAGTGTGCGCTGTGATGAGGCTGCTCTTTGAACCTGTCACTACACCTGTGGACATACGTTTTGCC
 TCAATGCCTTGAGCGTTGCCTGGACCACGCCCTCACTGTCCACTGTGCAAGGACAAGCTTTCTGAGCT
 GTTGGCAACCAGAAATTTAATGTAAGTGTCTGACTGAAGAGCTGATATCCGATACTTACCAGATGAG
 TTATCTGACAGGAAGAGGCTATGATGAAGAGATGTCAGAGCTGTCCAATCTGACCAGAGATGTGCCCA
 TCTTTGTGTGCGCATGGCCTTCCCTACGGTTCTTGCCTCTCCATGTGTTTCGAGCCCGCTACCGGCT
 CATGATAAGGAGATGCATGAAACTGGCACCAGCGGTTTGGCATGTGTTATCCGCGGAGAACGAGGG
 ATTTTCAGAGTATGGGTGCATGTTGGAGATAAAGGATGTCAGGACTTTCCCTGATGGAAGCTCTGTTGTGG
 ATGCAATTGGGATCAGCCGCTTCAGAGTCTGAGCCACCGCCACCGAGATGGCTACAACACGGCGGACAT
 CGAGTACCTGGAGGACGAGAAGGTGGAAGGTCCAGAGTTTGAAGAACTTACCGCCCTCCATGAGTCTGTT
 TACCAGCAGTCTGTGCTCTGGTTTGTCTCACTCCAGGATCACATGAAGAAACAAATCTTAGCCATTTTG
 GGTGATGCCAGACAGGAGCCAGGCCACAGACAATTCTAGTGGCCCTGCCTGGTCTGGTGGATTCT
 GGCCGTGTTGCCCTTGAACGCAAGGCTCAACTGGCATTCTCGGCATGGCCTCCCTGAAAGAGCGTCTT
 CTGGCCATTGACGAATATTAGTCATCATTACACGTAAGCTGAATAGTCGCCAAGAAATGGCTAACAAACA
 CACAGAGAGATAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR221330 representing NM_001029878
 Red=Cloning site Green=Tags(s)


MKNYDQALQDAEAVCQREPLLTKGHHIKAQILSGLGRHREVLKEFIYCLALNPECNSAKKETQKVICEVF
 FSTSESEHQTSTSTETGPEALCEGQTNQYPLEEAGGHANADNPKTPSEKSDAPADINSSVLYFILGLH
 CEEDKKALEGIVPAAPSSTLKRQLPSDAQDDEELKANTPEKIPKKDADSPQNRNASSLEEPEFTIDATD
 FECALCMRLLFEPVTPCGHTFCLKCLERCLDHAPHCLCKDKLSELLATRNFNVTVLTEELIFRYLPDE
 LSDRKRVYDEEMSELNLTRDVPIFVCAMAFPTVPCPLHVFEPYRLMIRRCMETGTRKFGMCLSAENAG
 ISEYGCMLKIDVTRFPDGSSVVDAGISRFRLSHRHRDGYNTADIEYLEDEKVEGPEFEELTALHESV
 YQSVSWFASLQDHMKQILSHFGSMPDREPEPQSNSSGPAWSWWILAVLPLERKAQLAILGMASLKERL
 LAIRRLVIITRKLNSRQEMANNTQRDN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9086_b03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001029878

ORF Size: 1554 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001029878.1](#), [NP_001025049.1](#)

RefSeq Size: 5077 bp

RefSeq ORF: 1557 bp
Locus ID: 381338
Cytogenetics: 1 B
MW: 58.4 kDa

Product images:

Circular map for MR221330