

Product datasheet for **MR221476L3V**

Adamts8 (NM_013906) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adamts8 (NM_013906) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Adamts8
Synonyms:	METH-2; METH2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013906
ORF Size:	2715 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221476).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013906.2 , NP_038934.2
RefSeq Size:	3632 bp
RefSeq ORF:	2718 bp
Locus ID:	30806
Cytogenetics:	9 16.49 cM



[View online »](#)

Gene Summary:

This gene encodes a member of "a disintegrin and metalloproteinase with thrombospondin motifs" (ADAMTS) family of multi-domain matrix-associated metalloendopeptidases that have diverse roles in tissue morphogenesis and pathophysiological remodeling, in inflammation and in vascular biology. This gene is expressed in mouse lung, heart and macrophage-rich areas of atherosclerotic plaques. The encoded preproprotein undergoes proteolytic processing to generate an active, zinc-dependent aggrecanase enzyme. This gene is located adjacent to a related ADAMTS gene on chromosome 9. [provided by RefSeq, May 2016]