

Product datasheet for **MR221520L4V**

Rab11fip1 (NM_029423) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rab11fip1 (NM_029423) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rab11fip1
Synonyms:	2010200K21Rik; 4833414G05Rik; Rcp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_029423
ORF Size:	1938 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221520).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029423.2 , NP_083699.2
RefSeq Size:	6402 bp
RefSeq ORF:	1938 bp
Locus ID:	75767
UniProt ID:	Q9D620
Cytogenetics:	8 A2
Gene Summary:	A Rab11 effector protein involved in the endosomal recycling process. Also involved in controlling membrane trafficking along the phagocytic pathway and phagocytosis (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)