

Product datasheet for MR221522

Klk15 (NM_174865) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Klk15 (NM_174865) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Klk15
Synonyms:	MGC129160; MGC129345
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR221522 representing NM_174865 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTGGCTTTCCTTGTTCCTTGTCTTGGTGTCTGCAGCGCAGGATGGCGACAAGGTGCTAGAAGGCG
AGGAGTGTGTGCCCCATTACAGCCCTGGCAAGTAGCCCTCTTCGAACGTGGCCGGTTAATTGTGGCGC
TTTCCTCATCTCACCACGCTGGGTGTGACTGCCGCCACTGCCAAACCCGTTTCATGAGAGTGCCTG
GGTGAACAATCTCCGCAAGTTCGATGGTCCAGAGCAACTGCGCTCCGTCTCTCGCATCATTCCACATC
CGGGCTACGAGGCTCGCACCCACCGCCATGACATTATGCTGTTGGGACTTTTCAAACCTGCCGGCTGAC
AGCCTACGTGCGTCTGTGGCGCTGCCAGACGTTGCCCTCTTATCGGTGAGGACTGCGTGGTGTGAGGC
TGGGGTCTGTTATCAGACAACAACCTGGAGCCACAGGGAGCCAGAAGTACACGTGAGACTCCCTGATA
CGTTGCATTGTGCCAACATCAGCATTATCTCTGAAGCATCTTGCAACAAGGACTACCCGGGTCGAGTGT
GCCACCATGGTGTGCGCGGCTCGAAGGCGGTGGCAGGACTCCTGTGAGGGTACTCCGGAGGACCC
TTGGTCTGTGGGGTGCCTGCAGGGCATTGTCTCTGGGGCGATGTCCCTGTGATACTACCAACG
CTGGCGTCTACAAAAAGTCTGCAGCTACTTGGAGTGGATCTGGGAAAATGTGAGAAGGTAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR221522 representing NM_174865
Red=Cloning site Green=Tags(s)


MWLLAFVLLVSAAQDGDKVLEGECEVPHSQPWQVALFERGRFNCGAFLISPRWVLTAAHCQTRFMRVRL
 GEHNLRFKFDGPEQLRSVSRIPHPGYEARTHHRDIMLLRFLFKPARLTAYVRPVALPRRCLIGEDCVVSG
 WGLLSDNPNPGATGSQKSHVRLPDTLHCANISIISEASCNKDYPGRVLPMTVCAGVEGGTDSCEGDSGGP
 LVCGGALQGI VSWGDVPCDTTTKPGVYTKVCSYLEWIWENVRRY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9007_b08.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_174865

ORF Size: 762 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_174865.2](#)

RefSeq Size: 765 bp


RefSeq ORF: 765 bp

Locus ID: 317652

Cytogenetics: 7 28.41 cM

MW: 28.5 kDa

Product images:



Circular map for MR221522