

## Product datasheet for MR221530L4

### Kctd17 (NM\_001081367) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Kctd17 (NM_001081367) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Kctd17
Synonyms:	2900008M13Rik; AA414907; AW742389; N28155
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221530).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



Locus ID: 72844

Cytogenetics: 15 E1