

## Product datasheet for **MR221561L4V**

### **Pax4 (NM\_001159926) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Pax4 (NM_001159926) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pax4
Synonyms:	Pax-4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001159926
ORF Size:	1011 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR221561).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001159926.1</a> , <a href="#">NP_001153398.1</a>
RefSeq Size:	1336 bp
RefSeq ORF:	1008 bp
Locus ID:	18506
UniProt ID:	<a href="#">P32115</a>
Cytogenetics:	6 11.99 cM
Gene Summary:	Plays an important role in the differentiation and development of pancreatic islet beta cells. Transcriptional repressor that competes with PAX6 in binding to a common element in the glucagon, insulin and somatostatin promoters.[UniProtKB/Swiss-Prot Function]



[View online »](#)